

zippie®
Explore your world!

Simba™

Generation 2015



A lightweight, energy efficient wheelchair with growth and adjustability





Simba™

Generation 2015




Innovative and adjustable, the wheelchair with growth!

One thing we're certain of is that kids will grow! That's why the Zippie Simba features 4cm of depth growth built into the frame, and 2cm of width growth without additional parts.

The Simba – your loyal friend through the years!



6 upholstery colours



A DESIGN THAT'S FIT FOR KIDS

With independence comes confidence, so every part of this wheelchair has been designed to help children be as independent as possible.

The Simba features a host of options geared towards energy efficiency and propulsion. With ample camber and a centre of gravity setting, children can experience a wheelchair that is easy to push. The abducted frame and safari brakes ensure that the child remains stable, even whilst interacting with their environment and playing with friends.



A FULLY ADJUSTABLE AXLE ADAPTER

The optimised axle adapter offers full adjustability of the rear seat height, centre of gravity and camber without compromising on any performance.



SIMBA - AN INNOVATIVE DESIGN

The outstanding design of the new castor link tube, combined with the innovative carbotechnology forks that are made from extremely lightweight stiff high-tech material, makes the wheelchair just that little bit lighter and easier to push!



THE CONVINCING SAFARI WHEEL LOCK

With its improved handling the safari side guard and wheel lock let you stop quickly and safely.



IT'S ALL ABOUT ERGONOMICS

The push bar can be easily adjusted in angle and height for an enhanced ergonomic position. The new stabilizer reinforces the backrest.

COLOURFUL OPTIONS

We want the Simba to be an expression of your child's unique personality. Aside from its super stylish frame, the Simba also offers 28 bright frame colours, six different upholstery colours, coloured castor forks and a full range of stylish spoke guards.

28 frame colours



TECHNICAL DATA

| | |
|-----------------------------|--|
| Maximum user weight: | 65kg |
| Seat width: | 22 – 36cm (growth +2cm without additional parts) |
| Seat depth: | 24 – 38cm (growth +4cm without additional parts) |
| Front seat height: | 36 – 48cm |
| Rear seat height: | 29 – 46cm |
| Backrest height: | 20 – 40cm |
| Backrest angle: | -16° to +20° |
| Frame versions: | Straight or abduction frame (3cm each side) |
| Frame angle: | 85° |
| Front wheel | 3"/ 4"/ 5"/ 6" |
| Rear wheel | 20"/ 22"/ 24" |
| Wheelchair weight: | From 8kg |
| Colours: | 28 colours |
| Crash tested: | Yes |



For further information on product options please check the order form on our website. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. You can view this document also in a larger font on a PDF file at www.SunriseMedical.co.uk



Sunrise Medical Limited

Thorns Road
 Brierley Hill
 West Midlands
 DY5 2LD
 Tel.: +44 (0) 845 605 66 88
 Fax: +44 (0) 845 605 66 89
 E-mail: enquiries@sunmed.co.uk
www.SunriseMedical.co.uk