

## Skippi

The Skippi children's power wheelchair gives the child an element of independence, giving them the opportunity to explore the world around them with their family and friends.

As the Skippi is an easily manoeuvrable, the child can explore every corner with ease. The Skippi can also be disassembled in a few easy steps to fit in to the boot of a small car for transport.

### Standard Model

- VR2 control system
- 6km/h speed
- 28Ah (C5) AGM batteries
- Battery charger
- 8" PU caster wheels
- Standard seat with back bag
- Seat cushion
- Telescoping side panels
- Individual telescoping leg supports
- Wheel cover for drive wheels
- Drive wheels with aluminium rim, black
- Lap belt with buckle



## Skippi

### Options

- 20° power seat tilt
- Power back support angle adjustment 30°
- Power seat height adjustment 20cm
- Leg support with single-panel aluminium footplate
- Height-adjustable, removable push handles
- Attendant control
- Tray
- Removable bumper bar
- LED lighting
- Choice of 5 frame colours



Attendant controller optional



Seat tilt optional



Seat riser optional



Transportation brackets



Skippi disassembled

### Technical Data

Seat width	300 – 420mm
Seat height	445mm
Seat depth	290 – 370mm
Overall length	870 / 940mm with leg support
Overall width	590mm
Back support height	400 / 420 / 440 / 460mm
Back support angle	-9° to +30°
Arm support height	160 – 230mm
Seat tilt	20° mechanical or optional power
Seat inclination	0° / 3° / 6°
Lower leg length	150 – 380mm
Max charging time	12 hours
Driving distance approx.	25km
Obstacle negotiation capacity	50mm
Climbing ability	7° (12%)
Turning radius	680mm
Weight of chair	66.5kg
Max user weight	50kg

### Technical Data

Speed	6 km/h
Battery capacity	28 Ah (C5) AGM batteries
Range under test conditions approx.	25 km
Control unit	Skippi (enAble40) Skippi Plus (enable50)
Max user weight	50kg
Obstacle height the wheelchair is able to overcome (optional)	50mm
Climbing ability	12%
Length including footplates	850mm
Width	570mm
Turning Radius	520mm
Seat height (without cushion) at 3° seat tilt	430mm
Seat width	300 – 380mm
Seat depth	280 – 380mm
Back angle	-9 to + 30°
Seat tilt	25° mechanical / 20° powered
Back height	340 / 460mm
Armrest height	170 – 230mm
Lower leg length	150 – 380mm
Weight when empty	62kg
Charging time	12hr
Skippi Plus with electric seat height adjustment	200mm
Seat inclination	0° / 3° / 6°