



RETAIL

SOLUTION SPP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____


PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urthane liner.

Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ishcial well, mid-level contour.



Weight capacities
 79.37 kg: 10"-14" width
 158.75 kg: 15"-21" width,
 226.79 kg: 22" width+
 Overall weight: 2.4kg
 Overall height 3"
 10x10" up to 24x24"
 Odd sizes in 17" and 19"
 Custom sizes available

CUSHION

		ADD WIDTH	ADD DEPTH	
Solution SPP. 10"-20" W.		ST-SOLSPP	<input type="text"/>	£190
Solution SPP. 22"-24" W		ST-SOLSPP	<input type="text"/>	£210
Cushion Cover.	Coolcore.	SOLSPP-NBLK	<input type="text"/>	STD
	Mesh.	SOLSPP-MBLK	<input type="text"/>	FOC
	Stealth Tek.	SOLSPP-TBLK	<input type="text"/>	FOC
Fluid Proof Liner.	Urethane.	URL1709070	<input type="text"/>	STD
	Polyurethane.	PUL1709030	<input type="text"/>	£48
Extra Covers And Liners.	Coolcore Cover.	SOLSPP-NBLK	<input type="text"/>	£70
	Mesh Cover.	SOLSPP-MBLK	<input type="text"/>	£44
	Stealth Tek Cover.	SOLSPP-TBLK	<input type="text"/>	£50
	Urethane Liner.	URL1709070	<input type="text"/>	£40
	Polyurethane Liner.	PUL1709030	<input type="text"/>	£48

		Width										
		10	12	14	15	16	17	18	19	20	22	24
Depth	10											
	12											
	14											
	15											
	16											
	17											
	18											
	19											
	20											
	22											
	24											

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

10"W x 10"D. <i>DWR1010E188A.</i>	£34	16"W x 16"D. <i>DWR1010E191G.</i>	£34	19"W x 19"D. <i>MEC1804152-1919.</i>	£34
10"W x 12"D. <i>DWR1010E188C.</i>	£34	16"W x 18"D. <i>DWR1010E191I.</i>	£34	20"W x 16"D. <i>DWR1010E193G.</i>	£34
12"W x 10"D. <i>DWR1010E189A.</i>	£34	16"W x 20"D. <i>DWR1010E191K.</i>	£34	20"W x 18"D. <i>DWR1010E193L.</i>	£34
12"W x 12"D. <i>DWR1010E189C.</i>	£34	17"W x 17"D. <i>DWR1010E234H.</i>	£34	20"W x 20"D. <i>DWR1010E193K.</i>	£34
12"W x 14"D. <i>DWR1010E189E.</i>	£34	17"W x 19"D. <i>DWR1010E234J.</i>	£34	20"W x 22"D. <i>DWR1010E193M.</i>	£34
14"W x 12"D. <i>DWR1010E190C.</i>	£34	18"W x 16"D. <i>DWR1010E192G.</i>	£34	22"W x 18"D. <i>DWR1010E194I.</i>	£34
14"W x 14"D. <i>DWR1010E190E.</i>	£34	18"W x 18"D. <i>DWR1010E192I.</i>	£34	22"W x 20"D. <i>DWR1010E194K.</i>	£34
14"W x 16"D. <i>DWR1010E190G.</i>	£34	18"W x 20"D. <i>DWR1010E192K.</i>	£34	22"W x 22"D. <i>DWR1010E194M.</i>	£34
14"W x 18"D. <i>DWR1010E190I.</i>	£34	18"W x 22"D. <i>DWR1010E192M.</i>	£34	24"W x 20"D. <i>DWR1010E236K.</i>	£34
15"W x 18"D. <i>DWR1010E232H.</i>	£34	19"W x 17"D. <i>MEC1804152-1917.</i>	£34	24"W x 22"D. <i>DWR1010E236M.</i>	£34

22" width and greater cushions include a flat rigidizer as standard.

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

15"W x 18"D. <i>CUS142467.</i>	£56	18"W x 16"D. <i>CUS142466.</i>	£56	20"W x 18"D. <i>CUS142465.</i>	£56
16"W x 16"D. <i>CUS142467.</i>	£56	18"W x 18"D. <i>CUS142466.</i>	£56	20"W x 19"D. <i>CUS142465.</i>	£56
16"W x 17"D. <i>CUS142467.</i>	£56	18"W x 19"D. <i>CUS142466.</i>	£56	20"W x 20"D. <i>CUS142465.</i>	£56
16"W x 18"D. <i>CUS142467.</i>	£56	18"W x 20"D. <i>CUS142466.</i>	£56		

For use on manual wheelchairs only.

NOTES & SPECIAL REQUESTS

Special Order Solution SPP Cushion.

£350

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW



TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

SOLUTION SPP - 15/01/24 V1



V1: 14/06/2022