

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 1/10/2024

Spark - Transit Not Included _____ £2250
165 lb weight capacity Not available with upholstery

Frame

Seat Width

| | |
|----------------|-----|
| 10" [SPSWT010] | N/C |
| 11" [SPSWT011] | N/C |
| 12" [SPSWT012] | N/C |
| 13" [SPSWT013] | N/C |
| 14" [SPSWT014] | N/C |
| 15" [SPSWT015] | N/C |
| 16" [SPSWT016] | N/C |

Seat Depth

| | |
|----------------|-----|
| 12" [SPSDE012] | N/C |
| 13" [SPSDE013] | N/C |
| 14" [SPSDE014] | N/C |
| 15" [SPSDE015] | N/C |
| 16" [SPSDE016] | N/C |
| 17" [SPSDE017] | N/C |
| 18" [SPSDE018] | N/C |

Transit Option

| | |
|--|-----|
| Transit System [SPTRNOPT] | N/C |
| <i>No charge when selected with Spark - Transit Included. Not available with Removable Seat Pan.</i> | |
| Omit Transit [SPTRNOP2] | N/C |



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Seat Pan

| | |
|--|------|
| Removable Aluminum Seat Pan [SPACC004] | £242 |
| <i>Pan depth 1/2" shorter than Frame Depth</i> | |
| Omit [SPACC021] | N/C |

Backrest

Backrest

Back Posts are folding and angle adjustable in 4° increments.

| Back Post | Price | Height Range | Back Height | | | | |
|--|-------|--------------------|-------------|------|----|----|----|
| | | | Short | Tall | 13 | 14 | 15 |
| Stroller Handle Back Post [SPBOP032] | £73 | Short | 20 | | | | |
| | | Tall | 24 | | | | |
| Height Adjustable Straight with Push Handle [SPBOP033] | £73 | Short ¹ | 13 | 14 | 15 | 16 | 17 |
| | | Medium | 16 | 17 | 18 | 19 | 20 |
| | | Tall ² | 20 | 21 | 22 | 23 | 24 |

¹ Not available with Dynamic Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest

Backrest Options

| | |
|--|-------|
| Removable Stroller Handle Extension [SPBOP012] | £121 |
| <i>Not available with Height Adjustable Straight with Push Handle Back Posts.</i> | |
| Dynamic Back (E2398) [SPBOP014] | £1932 |
| <i>Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".</i> | |

Frame Type / Footrest

Frame Type

| | |
|-----------------------------------|-----|
| Swing Away [SPFRM001] | N/C |
| Super Low - Swing Away [SPFRM002] | £49 |

Footrest Hanger Type

| Hanger | 60° | 70° | 80° | 90° | | Retail |
|--------------------------------------|------------|------------|-------------------------|-------------------------|------------|--------|
| Extension Tube Mount | [SPHAN001] | | [SPHAN003] ¹ | | | N/C |
| Front Mount | | [SPHAN016] | [SPHAN017] ¹ | [SPHAN004] ¹ | | N/C |
| Pediatric PRO ELR | | | | | [SPHAN044] | £290 |
| Pro Elevating Leg Rests ² | | | | | [SPHAN045] | £290 |
| Omit Hangers | | | | | [SPHAN046] | N/C |

¹ Not available with 7" or 8" Casters

² Not available with Seat Widths less than 14"

PRO ELR Cover

Available only on PRO ELR.

Gel Pad Covers (Pair) [SPELRPAD] _____ £304

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

| Footrest | Extension | Seat to Footrest Length | | | | | | | | | | | | |
|--------------------------------------|---------------------|-------------------------|-------------------|-------------------|-----------------|-------------------|-------|-------|-------|-------|-----------------------|-------|-------|-------|
| | | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | | | | |
| Extension Tube Hangers (60°,70°,80°) | Short | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | | | | |
| | Medium | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | | | | |
| | Long | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | | | | |
| Front Mount Hangers ¹ | Standard Selections | | | | | | | | | | Requires 3" Extension | | | |
| | Short | 4" ⁵ | 4.5" ⁵ | 5" | 5.5" | 6" | 6.5" | | | | | 10" | 10.5" | 11" |
| | | 7" | 7.5" | 8" | 8.5" | 9" | 9.5" | | | | | 11.5" | 12" | 12.5" |
| | Medium | 4" ⁵ | 4.5" ⁵ | 5" | 5.5" | 6" | 6.5" | 7" | 7.5" | | | 12" | 12.5" | 13" |
| | | 8" | 8.5" | 9" | 9.5" | 10" | 10.5" | 11" | 11.5" | | | 13.5" | 14" | 14.5" |
| | Long | 5" | 5.5" | 6" | 6.5" | 7" | 7.5" | 8" | 8.5" | 9" | | 14" | 14.5" | 15" |
| 9.5" | | 10" | 10.5" | 11" | 11.5" | 12" | 12.5" | 13" | 13" | | 15.5" | 16" | 16.5" | |
| Peds PRO ELR ² | Low Pivot | 7.5" ³ | 8" ³ | 8.5" ³ | 9" ⁴ | 9.5" ⁴ | 10" | 10.5" | 11" | 11.5" | 12" | 12.5" | 13" | |
| | High Pivot | 7.5" ³ | 8" ³ | 8.5" ³ | 9" ⁴ | 9.5" ⁴ | 10" | 10.5" | 11" | 11.5" | | | | |
| PRO ELR Low Pivot Setting | Short | 12" - 17.5" | | | | | | | | | | | | |
| | Long | 15" - 21" | | | | | | | | | | | | |
| PRO ELR High Pivot Setting | Short | 10.5" - 15.5" | | | | | | | | | | | | |
| | Long | 13" - 19" | | | | | | | | | | | | |

¹ See Minimum Front Seat Height chart above.

² Comes standard with 3" Extension

³ Not available on chairs 13" wide or more

⁴ Not available on chairs 15" wide or more

⁵ Available with 90° only

3" Extension

Available only with Front Mount Hangers.

3" Extension [SPEXT008] _____ £53

3" Extension (Provided as Parts) [SPEXT018] _____ £49

| Minimum Front Seat Height with Front Mount Hangers | | | | | | |
|--|-------|--------|--------|--------|-------|-----|
| Hanger Length | Short | | Medium | | Long | |
| | 70° | 80/90° | 70° | 80/90° | 70° | 80° |
| Minimum Front Seat Height with Standard Selection | 13" | 13.5" | 14" | 15" | 16" | 17" |
| Minimum Front Seat Height with 3" Extension | 14.5" | 16" | 16.5" | 18" | 18.5" | 20" |

Footplates

| | |
|--|------|
| Composite [SPFPL001] | N/C |
| <i>Not available with Front Mount Hanger or Frame Widths less than 14".</i> | |
| Composite Angle Adjustable [SPFPL014] | £182 |
| <i>Not available with Front Mount Hanger or Frame Widths less than 15".</i> | |
| Aluminum Angle Adjustable [SPFPL002] | N/C |
| Aluminum Locking Multi-Angle Adjustable [SPFPL020] | £303 |
| <i>Available only with Front Mount Hanger.</i> | |
| One Piece Flip Up Angle Adjustable [SPFPL009] | £188 |
| <i>Not available with Front Mount Hanger, PRO ELR or Frame Widths less than 14".</i> | |

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate.

| | |
|---|-----|
| Heel Loops [SPACC008] | £37 |
| <i>Not available with Frame Widths less than 14".</i> | |
| Adjustable Heel Loops [SPACC029] | £61 |
| Omit Heel Loops [SPACC011] | N/C |

Hanger Release

If no selection is made Classic Hanger Release will be used.

| | |
|--|-----|
| Classic [SPREL001] | N/C |
| 4-Way - In-Line [SPREL002] | N/C |
| <i>Not available with PRO ELR or Peds PRO ELR.</i> | |
| 4-Way - Out [SPREL003] | N/C |
| <i>Not available with PRO ELR or Peds PRO ELR.</i> | |

Calf Straps

| | |
|--------------------------------------|------|
| Velcro Adjustable [SPACC002] | £37 |
| Padded Velcro Adjustable [SPACC003] | £67 |
| Bodypoint Padded AeroMesh [SPACC010] | £103 |

Residual Limb Support

| | Limb Support | Price | Omit | Cover | | Pad Size | | | | | |
|-------|--------------------|-------|------------------------------|-------------------|---------------------|------------------|-------------------|-------------------|-------------------|--------------------|--------------------|
| | | | | Smooth [SPRLS006] | Air Mesh [SPRLS007] | 5 x 8 [SPRLS008] | 7 x 10 [SPRLS009] | 9 x 10 [SPRLS010] | 9 x 14 [SPRLS011] | 11 x 10 [SPRLS012] | 11 x 14 [SPRLS013] |
| Left | Left Limb Support | £170 | Omit Left Hanger [SPRLS004] | Smooth [SPRLS006] | Air Mesh [SPRLS007] | 5 x 8 [SPRLS008] | 7 x 10 [SPRLS009] | 9 x 10 [SPRLS010] | 9 x 14 [SPRLS011] | 11 x 10 [SPRLS012] | 11 x 14 [SPRLS013] |
| Right | Right Limb Support | £170 | Omit Right Hanger [SPRLS005] | Smooth [SPRLS020] | Air Mesh [SPRLS021] | 5 x 8 [SPRLS014] | 7 x 10 [SPRLS015] | 9 x 10 [SPRLS016] | 9 x 14 [SPRLS017] | 11 x 10 [SPRLS018] | 11 x 14 [SPRLS019] |

Front Frame Height and Casters

Front Seat Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

| Caster Size | Front Seat Height | | | | | | | | | | | | | | |
|-------------|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-----|
| | Uses Super Low Front Frame ¹ | | | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" |
| 3" | 13.5" | 14" | | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | | | |
| 4" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | | 18.5" | 19" | 19.5" | 20" | |
| 5" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | |
| 6" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | |
| 7" | 17" | 17.5" | | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | |
| 8" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | | | |

¹ Price is 42€

Castors

| Caster Type | | Caster Size | | | | | |
|-------------|----------------------------|-------------|------------|------------|------------|------------|------------|
| | | Price | 4" | 5" | 6" | 7" | 8" |
| | .75" Rollerblade | N/C | [SPCAS004] | | | | |
| | .75" Lighted Rollerblade | £86 | [SPCAS005] | | | | |
| | 1" Poly | N/C | [SPCAS006] | [SPCAS010] | [SPCAS014] | [SPCAS019] | [SPCAS020] |
| | 1" Poly Aluminum | £61 | [SPCAS007] | [SPCAS011] | | | |
| | 1" Pneumatic | N/C | | | [SPCAS017] | | [SPCAS021] |
| | 1.5" Poly | N/C | [SPCAS032] | [SPCAS025] | [SPCAS026] | | |
| | 1.5" Soft Roll Aluminum | £61 | [SPCAS009] | [SPCAS013] | [SPCAS016] | | |
| | 2" Poly | N/C | | | [SPCAS024] | | |
| | 2" Pneumatic w/Foam Insert | £218 | | | | | [SPCAS023] |

Caster Housing

| | |
|--|-----|
| Standard Housing [SPCSH001] | N/C |
| Performance Housing [SPCSH002] | £37 |
| <i>Not available with 350 lb Weight Capacity</i> | |

| Maximum Seat Dump Reference (Front Seat Height minus Rear Seat Height) | | | | | | | | |
|---|------|------|------|-----|-----|------|------|------|
| Seat Depth | 12" | 13" | 14" | 15" | 16" | 17" | 18" | 19" |
| Maximum Dump Standard Caster Housing | 1.5" | 1.5" | 1.5" | 2" | 2" | 2" | 2.5" | 2.5" |
| Maximum Dump Performance Caster Housing | 2.5" | 2.5" | 3" | 3" | 3" | 3.5" | 3.5" | 4" |

Fork Stem

| | |
|--|------|
| Best Fit [SPCST006] | N/C |
| <i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i> | |
| Standard Fork Stem [SPCST001] | N/C |
| 3/4" Longer Fork Stem [SPCST002] | N/C |
| 1 1/2" Longer Fork Stem [SPCST003] | N/C |
| Frog Legs Suspension Fork [SPCST013] | £344 |
| <i>Available only with 3", 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.</i> | |

Fork Options

Frog Leg Options

| | |
|-------------------------------------|-----|
| User Weight Under 80 lbs [SPCST004] | N/C |
| User Weight 80-120 lbs [SPCST008] | N/C |
| User Weight 120-160 lbs [SPCST009] | N/C |

Caster Pin Locks

| | |
|---|------|
| Caster Pin Locks [SPCST005] | £194 |
| <i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i> | |

Wheel Locks

Wheel Lock

| | |
|-------------------------|-----|
| Push to Lock [SPWLK002] | N/C |
| Pull to Lock [SPWLK003] | N/C |
| Omit [SPWLK007] | N/C |

Wheel Lock Options

| | |
|---|------|
| Extension Handles (EO961) [SPWLK004] | £49 |
| <i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i> | |
| Grade Aids [SPWLK009] | £145 |
| <i>Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires</i> | |

Rear Frame Height and Wheels

Rear Seat Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Frame Depths not available with 24" wheel size.

| Wheel Size | Tire Profile | Rear Seat Height | | | | | | | | | | |
|------------|--------------|--------------------|--------------------|--------------------|--------------------|--------------------|-------|--------------------|-------|--------------------|--------------------|--------------------|
| | | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" |
| 18" | Low | 11" ² | 11.5" ¹ | 12" ² | 12.5" ¹ | 13" ² | 13.5" | 14" ² | 14.5" | 15" ² | 15.5" ² | 16" ² |
| | Full | 11.5" ² | 12" ¹ | 12.5" ² | 13" ¹ | ² | 14" | 14.5" ² | 15" | 15.5" ² | 16" ² | 16.5" ² |
| 20" | Low | 11.5" ² | 12" ¹ | 12.5" ² | 13" ¹ | | 14" | 14.5" ² | 15" | 15.5" ² | 16" ² | 16.5" ² |
| | Full | 12" ² | 12.5" ¹ | 13" ² | 13.5" ¹ | 14" ² | 14.5" | 15" ² | 15.5" | 16" ² | 16.5" ² | 17" ² |
| 22" | Low | ² | 13" ¹ | 13.5" ² | 14" ¹ | 14.5" ² | 15" | 15.5" ² | 16" | 16.5" ² | 17" ² | 17.5" ² |
| | Full | 13" ² | 13.5" ¹ | 14" ² | 14.5" ¹ | 15" ² | 15.5" | 16" ² | 16.5" | 17" ² | 17.5" ² | 18" ² |
| 24" | Low | 13.5" ² | 14" ¹ | 14.5" ² | 15" ¹ | 15.5" ² | 16" | 16.5" ² | 17" | 17.5" ² | 18" ² | 18.5" ² |
| | Full | 14" ² | 14.5" ¹ | 15" ² | 15.5" ¹ | 16" ² | 16.5" | 17" ² | 17.5" | 18" ² | 18.5" ² | 19" ² |

¹ Amputee Axle Plate not available with lowest two Rear Frame Height positions
² Requires Universal Axle Plate

One Arm Drive

| | |
|--|------|
| One Arm Drive - Quick Release [SPDRV007] | £846 |
| <i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.</i> | |

Axle Plates

No selection needed with One Arm Drive.

| | |
|--|------|
| Universal [SPAXP009] | N/C |
| Amputee [SPAXP003] | £170 |
| <i>Not available with Wheelchair Drive Ready Option.</i> | |

Drive Handrim Side

| | |
|-----------------------------|-----|
| Left Side Drive [SPDRV002] | N/C |
| Right Side Drive [SPDRV003] | N/C |

Drive Handrim Type

| | |
|--------------------------------------|------|
| Aluminum Drive Handrim [SPDRV004] | N/C |
| Poly Coated Drive Handrim [SPDRV005] | £194 |

Wheels

| | | Rear Wheel Size | | | | |
|-----------------|--------------------------|-----------------|------------|-------------|-------------|------------|
| | | Price | 18" | 20" | 22" | 24" |
| Rear Wheel Type | MAXX Spoke | N/C | [SPWH1118] | [SPWH1120] | [SPWH1122] | [SPWH1124] |
| | Mag (Black) [SPMWC005] | £121 | | [SPWHE1220] | [SPWHE1222] | [SPWHE124] |
| | Colored Mag ¹ | £182 | | [SPWHE1320] | [SPWHE1322] | |
| | Spinergy Spox | £544 | | | [SPWHE422] | [SPWHE424] |

¹ Mag Color Selection on Page 7

Spinergy Option

Select spoke color and handrim tab for Spinergy wheels.

| | |
|-------------------------|------------------------|
| Black Spokes [SPSSC001] | N/C |
| Red Spokes [SPSSC002] | Blue Spokes [SPSSC005] |
| White Spokes [SPSSC003] | £91 |

Handrims

No selection needed with One Arm Drive.

Standard Handrims

| | |
|--|------|
| Aluminum Anodized [SPHRM001] | N/C |
| Plastic Coated [SPHRM002] | £121 |
| Projection (8 Vertical Projections) [SPHRM003] | £285 |
| <i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i> | |
| Omit Handrim [SPHRM012] | N/C |

Ergonomic Handrim

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

| | |
|---------------------------------------|------|
| Ergo Lite Regular Grip [SPHRM005] | £272 |
| Ergo Standard Regular Grip [SPHRM008] | £272 |
| Ergo Lite Extra Grip [SPHRM006] | £272 |
| Ergo Standard Extra Grip [SPHRM009] | £272 |
| Ergo Lite Omit Grip [SPHRM004] | £272 |
| Ergo Standard Omit Grip [SPHRM007] | £272 |

Handrim Options

Micro-Grip

| | |
|--|-----|
| 3M Micro Grip [SPACC006] | £91 |
| <i>Not available on Plastic Coated or Projection Handrims.</i> | |

Handrim Mount

| | |
|---|-----|
| Close Mount [SPTAB001] | N/C |
| <i>If no selection is made, Standard Mount will be used</i> | |

Rear Tire Type

Full Profile

| | |
|---|------|
| Street Pneumatic 1 3/8" [SPTIR001] | N/C |
| Full Poly 1 3/8" [SPTIR002] | N/C |
| Pneumatic with Airless Insert 1 3/8" [SPTIR004] | £110 |

Low Profile

| | |
|---|------|
| Low Poly 1" [SPTIR003] | N/C |
| Sentinel Solid Tire 1" [SPTIR005] | £110 |
| <i>Not available on 18" wheel size or 20" and 22" Mag wheels.</i> | |
| High Pressure 1" [SPTIR006] | N/C |
| <i>Not available on Mag wheel or 18" wheel size.</i> | |

Axles

Axles

No selection needed with One Arm Drive.

| | |
|---|------|
| Quick Release [SPAXO003] | N/C |
| Quad Quick Release [SPAXO001] | £127 |
| <i>Not available with Mag or Spinergy wheels.</i> | |

Axle Position

Rear Wheel Camber

| | | | |
|-------------|-------------|-------------|-----|
| 0° [SPWC00] | 2° [SPWC02] | 4° [SPWC04] | N/C |
|-------------|-------------|-------------|-----|

Center of Gravity

No selection needed with Amputee Axle Plate

| | | | |
|-----------------|-------------------|--------------------|-----|
| Rear [SPCOG003] | Center [SPCOG002] | Forward [SPCOG001] | N/C |
|-----------------|-------------------|--------------------|-----|

Armrests

| Armrest Type | Length | | Price |
|---|--------|------------|-------|
| Height Adjustable T-Arm Tall (11.5"-15") | Desk | [SPARM039] | N/C |
| | Full | [SPARM040] | |
| Height Adjustable T-Arm (8.5"-13.5") ¹ | Desk | [SPARM003] | £257 |
| | Full | [SPARM004] | |
| Height Adjustable T-Arm Low (6.5"-10.5") ² | Desk | [SPARM022] | £257 |
| | Full | [SPARM023] | |

| Armrest Type | Length | | Price |
|--|--------|------------|-------|
| Angle Adjustable Locking Extendable Flip Up ⁵ | Desk | [SPARM010] | £257 |
| | Full | [SPARM011] | |
| Pediatric Height Adjustable T-Arm Desk (6.5"-10") | | [SPARM014] | £189 |
| Tubular Flip Up ³ | | [SPARM015] | N/C |
| Omit Armrests | | [SPARM007] | N/C |

¹ Not available with Frame Depths less than 14".

² Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.

³ Not available with Arm Pads.

Arm Pads

| | |
|--|-----|
| Standard Arm Pad [SPARM026] | N/C |
| Waterfall Arm Pad [SPARM020] | £25 |
| <i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i> | |
| Foam Grip [SPARM025] | N/C |
| <i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i> | |

Accessories


Removable Side Guards

Not available with Height Adjustable T-Arm Armrests.

| | |
|---|-----|
| Adult Composite Side Guards (14" - 16") [SPSDG002] | N/C |
| Pediatric Composite Side Guard (10" - 14") [SPSDG001] | N/C |

Pelvic Positioning Belts

| | |
|--|------|
| 1" Auto Buckle Pelvic Positioning Belt [SPPOB004] | £68 |
| 1" Padded Auto Buckle Pelvic Positioning Belt [SPPOB005] | £151 |
| 1 1/2" Auto Buckle Pelvic Positioning Belt [SPPOB001] | £68 |
| 1 1/2" Padded Auto Buckle Pelvic Positioning Belt [SPPOB002] | £151 |
| 2" Airline Buckle Pelvic Positioning Belt [SPPOB003] | £110 |
| Small Bodypoint Evoflex (16-24") [SPPOB007] | £188 |
| Medium Bodypoint Evoflex (19-29") [SPPOB008] | £188 |
| Large Bodypoint Evoflex (25-37") [SPPOB009] | £188 |
| Belt Mounting Band Clamp [SPACC054] | £103 |



Anti-tips

| | |
|---------------------------|------|
| Rear Anti-Tips [SPATP001] | £212 |
| Omit [SPATP000] | N/C |

Spoke Guards

| | |
|--|------|
| Clear Spoke Guards [SPSPG001] | £151 |
| <i>Not available with Mag wheel.</i> | |
| Black Spoke Guards [SPSPG002] | £151 |
| <i>Not available with 18" or 20" wheel sizes or Mag wheel.</i> | |
| Blue Monster Spoke Guards (1) [SPSPG003] | £182 |
| Robo Fun Spoke Guards (1) [SPSPG004] | £182 |
| Tutu Cute Spoke Guards (1) [SPSPG005] | £182 |
| Tuxedo Time Spoke Guards (1) [SPSPG008] | £182 |

¹ Available only with 20" and 22" wheel sizes. Not available with Mag wheel.

Other Accessories

| | |
|--|------|
| Tool Kit [SPACC005] | £20 |
| Removable Underseat Pouch [SPACC015] | £49 |
| <i>Not available with 10" - 13" Seat Widths.</i> | |
| Pediatric Ki Mobility Backpack [SPACC026] | £100 |
| Ki Mobility Backpack [SPACC009] | £100 |
| Black Neoprene Impact Guards [SPIMG001] | £86 |
| Cane and Crutch Holder [SPACC007] | £91 |
| Canopy with Hardware [SPACC042] | £272 |
| <i>See image to the left</i> | |

Frame Color

Gloss Finish

| | |
|----------------------------|-----|
| Black [SPFCS001] | N/C |
| Blue Skies [SPFCS018] | N/C |
| Bubblegum Pink [SPFCS011] | N/C |
| Candy Blue [SPFCS012] | N/C |
| Candy Purple [SPFCS010] | N/C |
| Candy Red [SPFCS003] | N/C |
| Copper [SPFCS009] | N/C |
| Granny Smith [SPFCS008] | N/C |
| Metallic Orange [SPFCS032] | N/C |
| Yellow [SPFCS029] | N/C |

Matte Finish

| | |
|--------------------------------|-----|
| Military Green [SPFCS025] | N/C |
| Charcoal Gray [SPFCS031] | N/C |
| Emerald Green [SPFCS019] | N/C |
| Iridescent Burgundy [SPFCS021] | N/C |
| Iridescent Navy [SPFCS016] | N/C |
| Retro Red [SPFCS022] | N/C |
| Royal Blue [SPFCS023] | N/C |
| Shadow Black [SPFCS015] | N/C |
| Smooth White [SPFCS024] | N/C |

Hydrographics

| | |
|---|------|
| AOR1 (Desert Digital Camo) [SPGIT002] | £725 |
| Pink Camo [SPGIT005] | £725 |
| Muddy Girl [SPGIT010] | £725 |
| NEXT Camo [SPGIT003] | £725 |
| Carbon Fiber [SPGIT008] | £725 |
| <i>Select Hydrographics Base Color.</i> | |
| Zebra [SPGIT007] | £725 |
| <i>Select Hydrographics Base Color. Not available in Black.</i> | |

Hydrographics Base Color

| | | | |
|------------------|-------------------|-------------------------|-----|
| Black [SPGBC001] | Blue [SPGBC002] | Granny Smith [SPGBC004] | |
| Pink [SPGBC005] | Purple [SPGBC006] | Red [SPGBC007] | N/C |

Mag Wheel Color

Gloss Finish

| | | | | |
|----------------------------|-----------------------|---------------------------|-------------------------|-----|
| Black [SPMWC006] | Blue Skies [SPMWC007] | Bubblegum Pink [SPMWC008] | Candy Blue [SPMWC009] | |
| Candy Purple [SPMWC010] | Candy Red [SPMWC011] | Copper [SPMWC012] | Granny Smith [SPMWC013] | |
| Metallic Orange [SPMWC014] | Yellow [SPMWC017] | | | N/C |

Matte Finish

| | | | | |
|----------------------------|--------------------------|--------------------------|--------------------------------|-----|
| Military Green [SPMWC022] | Charcoal Gray [SPMWC018] | Emerald Green [SPMWC019] | Iridescent Burgundy [SPMWC020] | |
| Iridescent Navy [SPMWC021] | Retro Red [SPMWC023] | Royal Blue [SPMWC024] | Shadow Black [SPMWC025] | |
| Smooth White [SPMWC026] | | | | N/C |

| Spark Front Frame Height Matrix | | | | | | | | | | | | | | | | | | | | | | | |
|---------------------------------|-----------|-------|------|-------------------|-----------------|-------|------|-----------|-----------|-------|------|-----------|-----------|-------|-------|---|-----------|-------|------|-----------|-----------|-------|-------|
| 3" Caster | | | | 4" Caster | | | | 5" Caster | | | | 6" Caster | | | | 7" Caster | | | | 8" Caster | | | |
| FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork | FSH | Super Low | Stem | Fork |
| 14.5 | 13.5 | Std | 4" | 15 | 13.5 | Std | 4" | 16 | 14.5 | Std | 4" | 17 | 15.5 | Std | 5" | 18 | 16.5 | Std | 6" | 19 | 17.5 | Std | 6" |
| 15 | 14 | | | 15.5 | 14 | | | 16.5 | 15.5 | | | 3/4 | 17.5 | | | 16.5 | 3/4 | | | 18.5 | 17 | | |
| 15.5 | 14 | 3/4 | 4" | 15.5 ¹ | 14 ¹ | 3/4 | 4" | 17.5 | 16 | 1 1/2 | 5" | 18.5 | 17 | 1 1/2 | 6" | 18.5 | 17.5 | 3/4 | 6" | 20.5 | n/a | 1 1/2 | 7" |
| 16 | 14.5 | | | 16 | 14.5 | | | 16 | 14.5 | | | Std | 17 | | | 15.5 | Std | | | 17.5 | 16 | | |
| 16 | 15 | 1 1/2 | 4" | 16.5 ¹ | 15 ¹ | 1 1/2 | 5" | 16.5 | 15 | 3/4 | 5" | | 18 | 16.5 | 3/4 | 6" | | 19.5 | 18 | 1 1/2 | 7" | 20 | 18.5 |
| 16.5 | 15.5 | | | 17 | 15.5 | | | 16.5 | 15.5 | | | 3/4 | 17 | 16 | | | 1 1/2 | 6" | 17.5 | | | 16 | 3/4 |
| 15.5 | 14.5 | Std | 5" | 15.5 | 14 | Std | 5" | 17 | 16 | 3/4 | 6" | | 17.5 | 16 | 1 1/2 | 6" | | | 18 | 17 | 3/4 | 7" | |
| 16 | 15 | | | 16 | 14.5 | | | 16 | 15 | | | 3/4 | 16.5 | 15.5 | | | 1 1/2 | 6" | 18 | 16.5 | | | 1 1/2 |
| 16.5 | 15 | 3/4 | 5" | 16 | 15 | 3/4 | 5" | 16.5 | 15 | Std | 6" | | 18.5 | 17 | 1 1/2 | 7" | | | 19 | 18 | 3/4 | 7" | |
| 17 | 15.5 | | | 16.5 | 15.5 | | | 16.5 | 15 | | | 1 1/2 | 6" | 17 | | | 15.5 | Std | 6" | 18.5 | | | 17 |
| 17 | 16 | 1 1/2 | 5" | 17 | 15.5 | 1 1/2 | 6" | 17 | 15.5 | Std | 6" | | | 19 | 17.5 | 1 1/2 | 7" | | | 19 | 18 | 3/4 | 7" |
| 17.5 | 16.5 | | | 17.5 | 16 | | | 17.5 | 16 | | | 3/4 | 6" | 17.5 | 16 | | | 3/4 | 6" | 18 | 16.5 | | |
| 16.5 | 15.5 | Std | 6" | 16.5 | 15 | Std | 6" | 17.5 | 16 | 3/4 | 6" | | | 18 | 16.5 | 3/4 | 7" | | | 18.5 | 17 | Std | 20.5 |
| 17 | 16 | | | 17 | 15.5 | | | 17 | 16 | | | 3/4 | 6" | 18 | 16.5 | | | 1 1/2 | 6" | 18.5 | 17 | | 1 1/2 |
| 17.5 | 16 | 3/4 | 6" | 17 | 16 | 3/4 | 6" | 18.5 | 17 | 1 1/2 | 6" | | | 19 | 17.5 | 1 1/2 | 7" | | | 19.5 | 18 | 1 1/2 | |
| 18 | 16.5 | | | 17.5 | 16.5 | | | 17.5 | 16.5 | | | 1 1/2 | 6" | 18 | 16.5 | | | 1 1/2 | 6" | 18.5 | 17 | | 1 1/2 |
| 18 | 17 | 1 1/2 | 6" | 18 | 16.5 | 1 1/2 | 6" | 18.5 | 17 | 1 1/2 | 6" | | | 18.5 | 17.5 | 3/4 | 7" | | | 20 | 18.5 | 1 1/2 | |
| 18.5 | 17.5 | | | 18.5 | 17 | | | 18.5 | 17 | | | 1 1/2 | 6" | 19 | 17.5 | | | 1 1/2 | 6" | 19 | 18 | | 3/4 |
| | | 18 | 16.5 | 18 | 16.5 | Std | 7" | 18 | 16.5 | Std | 7" | | | 19.5 | 18.5 | 1 1/2 | 7" | | | 20 | 19 | 1 1/2 | |
| | | 18.5 | 17 | 18.5 | 17 | | | 3/4 | 7" | | | 18.5 | 17 | Std | 7" | | | 20 | 19 | 1 1/2 | 7" | | 20 |
| | | 18.5 | 17.5 | 19 | 17.5 | 1 1/2 | 7" | | | 19 | 17.5 | 1 1/2 | 7" | | | 19.5 | 18 | 1 1/2 | 7" | | | 20 | 18.5 |
| | | 19 | 18 | 18.5 | 17.5 | | | 1 1/2 | 7" | 18.5 | 17.5 | | | 3/4 | 7" | 20 | 18.5 | | | 1 1/2 | 7" | 20 | 18.5 |
| | | 19.5 | 18.5 | 19 | 18 | 1 1/2 | 7" | | | 19 | 18 | 3/4 | 7" | | | 20.5 | 19 | 1 1/2 | 7" | | | 20.5 | 19 |
| | | 20 | 19 | 19.5 | 18.5 | | | 1 1/2 | 7" | 19.5 | 18 | | | 1 1/2 | 7" | 21 | n/a | | | 1 1/2 | 7" | 21 | n/a |
| | | | | | | | | | | 19.5 | 18 | 1 1/2 | 7" | | | 6" x 2" caster not available with 5" fork | | | | | | | |
| | | | | | | | | 20 | 18.5 | 1 1/2 | 7" | | | | | | | | | | | | |
| | | | | | | | | 20.5 | 19 | | | | | | | | | | | | | | |

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Axiom G

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover |
| Retail | | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom S

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover |
| Retail | | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom P

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover |
| Retail | | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom SP - Visco

| | | Width | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|--|-----|--|--|
| | | 14" | | | 15" | | | 16" | | | 17" | | | 18" | | | 19" | | | 20" | | | 21" | | | 22" | | |
| | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | | | |
| Retail | | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | | | | | |
| Height | 14" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Axiom SP - Fluid

| | | Width | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|--|-----|--|--|
| | | 14" | | | 15" | | | 16" | | | 17" | | | 18" | | | 19" | | | 20" | | | 21" | | | 22" | | |
| | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | | | |
| Retail | | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | | | | | |
| Height | 14" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Solid Seat Insert

Weight capacity: 350 lbs

| | | Width | | | | | | | | |
|--------|-----|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| | | 14" | 15" | 16" | 17" | 18" | 19" | 20" | 21" | 22" |
| | | Cushion |
| Retail | | £55 | £55 | £55 | £55 | £55 | £55 | £55 | £55 | £55 |
| Height | 14" | | | | | | | | | |
| | 15" | | | | | | | | | |
| | 16" | | | | | | | | | |
| | 17" | | | | | | | | | |
| | 18" | | | | | | | | | |
| | 19" | | | | | | | | | |
| | 20" | | | | | | | | | |
| | 21" | | | | | | | | | |
| | 22" | | | | | | | | | |

| Axiom Back Width | Fits Back Frame Widths | |
|------------------|-----------------------------|-----------------------------|
| | Inner Slots | Outer Slots |
| 12" (30.5cm) | 11" (28.0cm) | 12" -14" (30.5cm -35.5cm) |
| 14" (35.5cm) | 12" -13" (30.5cm -33.0cm) | 14" -16" (35.5cm -40.5cm) |
| 16" (40.5cm) | 14" -15" (35.5cm -38.0cm) | 16" -18" (40.5cm -45.5cm) |
| 18" (45.5cm) | 16" - 17" (40.5cm - 43.0cm) | 18" - 20" (45.5cm - 50.5cm) |
| 20" (51.0cm) | 18" -19" (46.0cm -48.5cm) | 20" -22" (50.5cm -55.5cm) |
| 22" (56.0cm) | 20" -21" (51.0cm -53.5cm) | 22" -24" (56.0cm -61.0cm) |

| Hardware Description | Code |
|--|------------|
| Improved DEX Hardware, removable adjustable hardware (300 lb max user weight) | DEX |
| New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight) | BTX |
| New Fixed Hardware, non-removeable hardware up to (300 lb max user weight) | FXT |



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

| | | Width | | | | | | |
|--------|-----|-------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| Height | Low | 8" | [AXPB1208] | [AXPB1408] | [AXPB1608] | [AXPB1808] | [AXPB2008] | |
| | | 10" | | [AXPB1410] | [AXPB1610] | [AXPB1810] | [AXPB2010] | [AXPB2210] |
| | | 12" | [AXPB1212] | [AXPB1412] | [AXPB1612] | [AXPB1812] | [AXPB2012] | |
| | Mid | 14" | | [AXPB1414] | [AXPB1614] | [AXPB1814] | [AXPB2014] | [AXPB2214] |
| | | 16" | [AXPB1216] | [AXPB1416] | [AXPB1616] | [AXPB1816] | [AXPB2016] | |
| | | 18" | | [AXPB1418] | [AXPB1618] | [AXPB1818] | [AXPB2018] | [AXPB2218] |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

| Lateral Trunk Supports for AXPB Backs | | | | | | | |
|---------------------------------------|------|-------|---|-----------|------|------|--|
| Pair | Left | Right | Description | Item Code | Pair | Each | |
| | | | 3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back | PBLF3x4P | £295 | £148 | |
| | | | 4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back | PBLF4x5P | £295 | £148 | |
| | | | 3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back | PBLS3x4P | £628 | £314 | |
| | | | 4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back | PBLS4x5P | £628 | £314 | |
| | | | 3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware | PBLAS3x4P | £628 | £314 | |
| | | | 4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware | PBLAS4x5P | £628 | £314 | |

Axiom Posterior Lateral Back (PL)

| | | Width | | | | | | |
|--------|-----|------------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| Height | Low | 8" | [AXPL1208] | [AXPL1408] | [AXPL1608] | [AXPL1808] | [AXPL2008] | |
| | | 10" | | [AXPL1410] | [AXPL1610] | [AXPL1810] | [AXPL2010] | [AXPL2210] |
| | | 12" | [AXPL1212] | [AXPL1412] | [AXPL1612] | [AXPL1812] | [AXPL2012] | |
| | Mid | 14" | | [AXPL1414] | [AXPL1614] | [AXPL1814] | [AXPL2014] | [AXPL2214] |
| | | 16" | [AXPL1216] | [AXPL1416] | [AXPL1616] | [AXPL1816] | [AXPL2016] | |
| | | 18" | | [AXPL1418] | [AXPL1618] | [AXPL1818] | [AXPL2018] | [AXPL2218] |
| Tall | 20" | [AXPL1220] | [AXPL1420] | [AXPL1620] | [AXPL1820] | [AXPL2020] | | |
| | 22" | | [AXPL1422] | [AXPL1622] | [AXPL1822] | [AXPL2022] | [AXPL2222] | |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

| Lateral Trunk Supports for AXPL & AXLB Backs | | | | | | |
|--|------|-------|---|-----------|------|------|
| Pair | Left | Right | Description | Item Code | Pair | Each |
| | | | 3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back | PLLF3x4P | £295 | £295 |
| | | | 4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back | PLLF4x5P | £295 | £295 |
| | | | 3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back | PLLS3x4P | £628 | £314 |
| | | | 4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back | PLLS4x5P | £628 | £314 |
| | | | 3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware | PLLAS3x4P | £628 | £314 |
| | | | 4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware | PLLAS4x5P | £628 | £314 |

Axiom Deep Lateral Back (DL)

| | | Width | | | | | | |
|--------|------|-------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| HCPC | | E2620 | E2620 | E2620 | E2620 | E2620 | E2621 | |
| Height | Mid | 14" | [AXDL1214] | [AXDL1414] | [AXDL1614] | [AXDL1814] | [AXDL2014] | [AXDL2214] |
| | | 16" | | [AXDL1416] | [AXDL1616] | [AXDL1816] | [AXDL2016] | |
| | | 18" | [AXDL1218] | [AXDL1418] | [AXDL1618] | [AXDL1818] | [AXDL2018] | [AXDL2218] |
| | Tall | 20" | [AXDL1220] | [AXDL1420] | [AXDL1620] | [AXDL1820] | [AXDL2020] | |
| | | 22" | | [AXDL1422] | [AXDL1622] | [AXDL1822] | [AXDL2022] | [AXDL2222] |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

| Description | Std. Straight Removable Adjustable Hardware | | Multi-Axis Removable Adjustable Hardware | | Posalinc Removable Adjustable Hardware | | HCPCS |
|-------------------------------------|---|--------|--|--------|--|--------|---------------|
| | Item Code | Retail | Item Code | Retail | Item Code | Retail | |
| Comfort Soft Headrest 8" | [CS08SS] | £268 | [CS08MA] | £268 | [CS08PL] | £268 | E0955 / E1028 |
| Comfort Soft Headrest 10" | [CS10SS] | £268 | [CS10MA] | £268 | [CS10PL] | £268 | E0955 / E1028 |
| Comfort Soft Headrest 14" | [CS14SS] | £268 | [CS14MA] | £268 | [CS14PL] | £268 | E0955 / E1028 |
| Contour Positioning Pad, 8" X 6.25" | [CP08SS] | £303 | [CP08MA] | £303 | [CP08PL] | £303 | E0955 / E1028 |
| Contour Positioning Pad, 9" X 6.5" | [CP09SS] | £303 | [CP09MA] | £303 | [CP09PL] | £303 | E0955 / E1028 |
| Contour Positioning Pad, 10" X 7.5" | [CP10SS] | £303 | [CP10MA] | £303 | [CP10PL] | £303 | E0955 / E1028 |
| Tri-soft Headrest 10" | [TS10SS] | £386 | [TS10MA] | £386 | [TS10PL] | £386 | E0955 / E1028 |
| Tri-soft Headrest 14" | [TS14SS] | £386 | [TS14MA] | £386 | [TS14PL] | £386 | E0955 / E1028 |

| Description | Item Code | Retail | HCPCS |
|-----------------------------------|-----------|--------|-------|
| Universal Headrest Mounting Plate | [UMP] | £43 | K0108 |

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports

| Description | Item Code | HCPCS | Retail |
|--|-----------|-------|--------|
| Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small | [14941] | E0960 | £100 |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium | [14942] | E0960 | £100 |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large | [14943] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small | [14924] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium | [14927] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large | [14928] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small | [14851] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Small | [14852] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Medium | [14853] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Large | [14854] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small | [14856] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small | [14857] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium | [14858] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large | [14859] | E0960 | £100 |
| AirLogic Anterior Trunk Support: Standard, Structure, Small | [14282] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, Medium | [14283] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, Large | [14284] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, X-Large | [14285] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small | [14257] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium | [14258] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large | [14259] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large | [14260] | E0960 | £300 |

Accessories

| | |
|--|-----|
| Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i> | £43 |
| Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i> | N/C |
| Privacy Flap (K0108) [PRFLAP] | £44 |

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

| | |
|---|-----|
| Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i> | N/C |
|---|-----|

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

| | |
|-------------------------------------|-----|
| Black Back Shell [AXMBLK] | N/C |
| Blue Skies Back Shell [AXMSKY] | N/C |
| Bubblegum Pink Back Shell [AXMPNK] | N/C |
| Candy Blue Back Shell [AXMCBL] | N/C |
| Candy Purple Back Shell [AXMPUR] | N/C |
| Candy Red Back Shell [AXMRED] | N/C |
| Copper Back Shell [AXMCOP] | N/C |
| Granny Smith Back Shell [AXMGSM] | N/C |
| Metallic Orange Back Shell [AXMMOR] | N/C |
| Yellow Back Shell [AXMYLW] | N/C |
| Silver Back Shell [AXMSIL] | N/C |
| Teal Back Shell [AXMTEA] | N/C |

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

| | |
|---|-----|
| Military Green Back Shell [AXMAGM] | N/C |
| Charcoal Gray Back Shell [AXMCGM] | N/C |
| Emerald Green Back Shell [AXMEGM] | N/C |
| Iridescent Burgundy Back Shell [AXMIBM] | N/C |
| Iridescent Navy Back Shell [AXMINM] | N/C |
| Retro Red Back Shell [AXMRRM] | N/C |
| Royal Blue Back Shell [AXMRBM] | N/C |
| Smooth White Back Shell [AXMSWM] | N/C |