

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 6/17/2024

Spark - Transit Not Included £2250

165 lb weight capacity Not available with upholstery

Frame

Seat Width

10" [SPSWT010]	N/C
11" [SPSWT011]	N/C
12" [SPSWT012]	N/C
13" [SPSWT013]	N/C
14" [SPSWT014]	N/C
15" [SPSWT015]	N/C
16" [SPSWT016]	N/C

Seat Depth

12" [SPSDE012]	N/C
13" [SPSDE013]	N/C
14" [SPSDE014]	N/C
15" [SPSDE015]	N/C
16" [SPSDE016]	N/C
17" [SPSDE017]	N/C
18" [SPSDE018]	N/C

Transit Option

Transit System [SPTRNOPT]	N/C
<i>No charge when selected with Spark - Transit Included</i>	
Omit Transit [SPTRNOP2]	N/C



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Seat Pan

Removable Aluminum Seat Pan [SPACC004]	£242
<i>Pan depth 1/2" shorter than Frame Depth</i>	
Omit [SPACC021]	N/C

Backrest

Backrest

Back Posts are folding and angle adjustable in 4° increments.

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post [SPBOP032]	£73	Short	20				
		Tall	24				
Height Adjustable Straight with Push Handle [SPBOP033]	£73	Short ¹	13	14	15	16	17
		Medium	16	17	18	19	20
		Tall ²	20	21	22	23	24

¹ Not available with Dynamic Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest

Backrest Options

Removable Stroller Handle Extension [SPBOP012]	£121
<i>Not available with Height Adjustable Straight with Push Handle Back Posts.</i>	
Dynamic Back (E2398) [SPBOP014]	£1932
<i>Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".</i>	

Frame Type / Footrest

Frame Type

Swing Away [SPFRM001]	N/C
Super Low - Swing Away [SPFRM002]	£49

Footrest Hanger Type

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	[SPHAN001]		[SPHAN003] ¹			N/C
Front Mount		[SPHAN016]	[SPHAN017] ¹	[SPHAN004] ¹		N/C
Pediatric PRO ELR ³					[SPHAN044]	£290
Pro Elevating Leg Rests ^{2,3}					[SPHAN045]	£290
Omit Hangers					[SPHAN046]	N/C

- ¹ Not available with 7" or 8" Casters
- ² Not available with Seat Widths less than 14"
- ³ Available only with Classic Hanger Release

Pro ELR Cover

Available only on PRO ELR.

Pro ELR Gel Pad Covers (PAIR) [SPELRPAD]	£304
------------------------------------------	------

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length											
		11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"			
Extension Tube Hangers (60°,70°,80°)	Short	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"			
	Medium	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"			
	Long	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"			
Front Mount Hangers ¹	Short Front Mount Hanger (70°, 80°, 90°)	Standard Selections									Requires 3" Extension		
		4" ⁴	4.5" ⁴	5"	5.5"	6"	6.5"				10"	10.5"	11"
	Medium Front Mount Hanger (70°, 80°, 90°)	7"	7.5"	8"	8.5"	9"	9.5"				11.5"	12"	12.5"
		4" ⁴	4.5" ⁴	5"	5.5"	6"	6.5"	7"	7.5"		12"	12.5"	13"
	Long Front Mount Hanger (70°, 80°)	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"		13.5"	14"	14.5"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14"	14.5"	15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"
Peds PRO ELR	Low Pivot	7.5"	8" ³	8.5" ³	9"	9.5"	10"	10.5" ²	11" ²	11.5" ²	12" ²	12.5" ²	13" ²
	High Pivot	7.5" ³	8" ³	8.5" ³	9" ²	9.5" ²	10" ²	10.5" ²	11" ²	11.5" ²			
PRO ELR Low Pivot Setting	Short	12" - 17.5"											
	Long	15" - 21"											
PRO ELR High Pivot Setting	Short	10.5" - 15.5"											
	Long	13" - 19"											

- ¹ See Minimum Front Seat Height chart above
- ² Comes standard with 3" Extension
- ³ Not available on chairs 14" wide or more
- ⁴ Available with 90° only

3" Extension

Available only with Front Mount Hangers.

3" Extension [SPEXT008]	£53
3" Extension (Provided as Parts) [SPEXT018]	£49

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Footplates

Composite [SPFPL001]	N/C
<i>Not available with Front Mount Hanger or Frame Widths less than 14".</i>	
Composite Angle Adjustable [SPFPL014]	£182
<i>Not available with Front Mount Hanger or Frame Widths less than 15".</i>	
Aluminum Angle Adjustable [SPFPL002]	N/C
Aluminum Locking Multi-Angle Adjustable [SPFPL020]	£303
<i>Available only with Front Mount Hanger.</i>	
One Piece Flip Up Angle Adjustable [SPFPL009]	£188
<i>Not available with Front Mount Hanger, PRO ELR or Frame Widths less than 14".</i>	

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate.

Heel Loops [SPACC008]	£37
<i>Not available with Frame Widths less than 14".</i>	
Adjustable Heel Loops [SPACC029]	£61
Omit Heel Loops [SPACC011]	N/C

Hanger Release

If no selection is made Classic Hanger Release will be used.

Classic [SPREL001]	N/C
4-Way - In-Line [SPREL002]	N/C
<i>Not available with PRO ELR or Peds PRO ELR.</i>	
4-Way - Out [SPREL003]	N/C
<i>Not available with PRO ELR or Peds PRO ELR.</i>	

Calf Straps

Velcro Adjustable [SPACC002]	£37
Padded Velcro Adjustable [SPACC003]	£67
Padded Airmesh Calf Strap [SPACC010]	£103

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
				Smooth [SPRLS006]	Air Mesh [SPRLS007]	5 x 8 [SPRLS008]	7 x 10 [SPRLS009]	9 x 10 [SPRLS010]	9 x 14 [SPRLS011]	11 x 10 [SPRLS012]	11 x 14 [SPRLS013]
Left	Left Limb Support	£170	Omit Left Hanger [SPRLS004]	Smooth [SPRLS006]	Air Mesh [SPRLS007]	5 x 8 [SPRLS008]	7 x 10 [SPRLS009]	9 x 10 [SPRLS010]	9 x 14 [SPRLS011]	11 x 10 [SPRLS012]	11 x 14 [SPRLS013]
Right	Right Limb Support	£170	Omit Right Hanger [SPRLS005]	Smooth [SPRLS020]	Air Mesh [SPRLS021]	5 x 8 [SPRLS014]	7 x 10 [SPRLS015]	9 x 10 [SPRLS016]	9 x 14 [SPRLS017]	11 x 10 [SPRLS018]	11 x 14 [SPRLS019]

Front Frame Height and Casters

Front Seat Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

Caster Size	Front Seat Height														
	Uses Super Low Front Frame ¹			14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
3"	13.5"	14"		14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"			
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"		18.5"	19"	19.5"	20"	
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"				
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"						
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								

¹ Price is 43€

Castors

Caster Type	Price	Caster Size				
		4"	5"	6"	7"	8"
.75" Rollerblade	N/C	[SPCAS004]				
.75" Lighted Rollerblade	£86	[SPCAS005]				
1" Poly	N/C	[SPCAS006]	[SPCAS010]	[SPCAS014]	[SPCAS019]	[SPCAS020]
1" Poly Aluminum	£61	[SPCAS007]	[SPCAS011]			
1" Pneumatic	N/C			[SPCAS017]		[SPCAS021]
1.5" Poly	N/C	[SPCAS032]	[SPCAS025]	[SPCAS026]		
1.5" Soft Roll Aluminum	£61	[SPCAS009]	[SPCAS013]	[SPCAS016]		
2" Poly	N/C			[SPCAS024]		
2" Pneumatic w/Foam Insert	£218					[SPCAS023]



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Caster Housing

Standard Housing [SPCSH001]	N/C
Performance Housing [SPCSH002]	£37
<i>Not available with 350 lb Weight Capacity</i>	

Maximum Seat Dump Reference (Front Seat Height minus Rear Seat Height)								
Seat Depth	12"	13"	14"	15"	16"	17"	18"	19"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"	2.5"
Maximum Dump Performance Caster Housing	2.5"	2.5"	3"	3"	3"	3.5"	3.5"	4"

Fork Stem

Best Fit [SPCST006]	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem [SPCST001]	N/C
3/4" Longer Fork Stem [SPCST002]	N/C
1 1/2" Longer Fork Stem [SPCST003]	N/C
Frog Legs Suspension Fork [SPCST013]	£344
<i>Available only with 3", 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.</i>	

Fork Options

Frog Leg Options

User Weight Under 80 lbs [SPCST004]	N/C
User Weight 80-120 lbs [SPCST008]	N/C
User Weight 120-160 lbs [SPCST009]	N/C

Caster Pin Locks

Caster Pin Locks [SPCST005]	£194
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Wheel Locks

Wheel Lock

Push to Lock [SPWLK002]	N/C
Pull to Lock [SPWLK003]	N/C
Omit [SPWLK007]	N/C

Wheel Lock Options

Extension Handles (EO961) [SPWLK004]	£49
Grade Aids [SPWLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires</i>	

Rear Frame Height and Wheels

Rear Seat Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Frame Depths not available with 24" wheel size.

Wheel Size	Tire Profile	Rear Seat Height										
		11" ²	11.5" ¹	12" ²	12.5" ¹	13" ²	13.5"	14" ²	14.5"	15"	15.5" ²	16" ²
18"	Low	11" ²	11.5" ¹	12" ²	12.5" ¹	13" ²	13.5"	14" ²	14.5"	15"	15.5" ²	16" ²
	Full	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5" ²	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
20"	Low	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5"	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
	Full	12" ²	12.5" ¹	13" ²	13.5" ¹	14" ²	14.5"	15" ²	15.5"	16" ²	16.5" ²	17" ²
22"	Low	12.5" ²	13" ¹	13.5" ²	14" ¹	14.5" ²	15"	15.5" ²	16"	16.5" ²	17" ²	17.5" ²
	Full	13" ²	13.5" ¹	14" ²	14.5" ¹	15" ²	15.5"	16" ²	16.5"	17" ²	17.5" ²	18" ²
24"	Low	13.5" ²	14" ¹	14.5" ²	15" ¹	15.5" ²	16"	16.5" ²	17"	17.5" ²	18" ²	18.5" ²
	Full	14" ²	14.5" ¹	15" ²	15.5" ¹	16" ²	16.5"	17" ²	17.5"	18" ²	18.5" ²	19" ²

¹ Amputee Axle Plate not available with lowest two Rear Frame Height positions
² Requires Universal Axle Plate

One Arm Drive

One Arm Drive - Quick Release [SPDRV007]	£846
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.</i>	

Axle Plates

No selection needed with One Arm Drive.

Universal [SPAXP009]	N/C
Amputee [SPAXP003]	£170
<i>Not available with Wheelchair Drive Ready Option.</i>	

Drive Handrim Side

Left Side Drive [SPDRV002]	N/C
Right Side Drive [SPDRV003]	N/C

Drive Handrim Type

Aluminum Drive Handrim [SPDRV004]	N/C
Poly Coated Drive Handrim [SPDRV005]	£194

Wheels

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel Type	MAXX Spoke	N/C	[SPWH1118]	[SPWH1120]	[SPWH1122]	[SPWH1124]
	Mag (Black) [SPMWC005]	£121		[SPWHE1220]	[SPWHE1222]	[SPWHE124]
	Colored Mag ¹	£182		[SPWHE1320]	[SPWHE1322]	
	Spinergy Spox	£544			[SPWHE422]	[SPWHE424]

¹ Mag Color Selection on Page 7

Spinergy Option

Select spoke color for Spinergy wheels.

Black Spokes [SPSSC001]	N/C
Red Spokes [SPSSC002]	Blue Spokes [SPSSC005]
White Spokes [SPSSC003]	£91

Handrims

Standard Handrims

Aluminum Anodized [SPHRM001]	N/C
Plastic Coated [SPHRM002]	£121
Projection (8 Vertical Projections) [SPHRM003]	£285
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels. Not available with One Arm Drive.</i>	
Omit Handrim [SPHRM012]	N/C

Ergonomic Handrim

Not available with 20" wheels

Ergo Lite Regular Grip [SPHRM005]	£272
Ergo Standard Regular Grip [SPHRM008]	£272
Ergo Lite Extra Grip [SPHRM006]	£272
Ergo Standard Extra Grip [SPHRM009]	£272
Ergo Lite Omit Grip [SPHRM004]	£272
Ergo Standard Omit Grip [SPHRM007]	£272

Handrim Options

3M Micro Grip [SPACC006]	£91
<i>Not available on Plastic Coated or Projection Handrims.</i>	

Handrim Mount

Close Mount [SPTAB001]	N/C
<i>If no selection is made, Standard Mount will be used. Reduces overall width by .75".</i>	

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [SPTIR001]	N/C
Full Poly 1 3/8" [SPTIR002]	N/C
Pneumatic with Airless Insert 1 3/8" [SPTIR004]	£110

Low Profile

Low Poly 1" [SPTIR003]	N/C
Sentinel Solid Tire 1" [SPTIR005]	£110
<i>Not available on 18" wheel size or 20" and 22" Mag wheels.</i>	
High Pressure 1" [SPTIR006]	N/C
<i>Not available on Mag wheel or 18" wheel size.</i>	

Axles

Axles

No selection needed with One Arm Drive.

Quick Release [SPAXO003]	N/C
Quad Quick Release [SPAXO001]	£127
<i>Not available with Mag or Spinergy wheels.</i>	

Axle Position

Rear Wheel Camber

0° [SPWC00]	2° [SPWC02]	4° [SPWC04]	N/C
-------------	-------------	-------------	-----

Center of Gravity

No selection needed with Amputee Axle Plate

Rear [SPCOG003]	Center [SPCOG002]	Forward [SPCOG001]	N/C
-----------------	-------------------	--------------------	-----

Armrests

Armrest Type	Length		Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk	[SPARM039]	N/C
	Full	[SPARM040]	
Height Adjustable T-Arm (8.5"-13.5") ¹	Desk	[SPARM003]	£257
	Full	[SPARM004]	
Height Adjustable T-Arm Low (6.5"-10.5") ²	Desk	[SPARM022]	£257
	Full	[SPARM023]	

Armrest Type	Length		Price
Angle Adjustable Locking Extendable Flip Up	Desk	[SPARM010]	£257
	Full	[SPARM011]	
Pediatric Height Adjustable T-Arm Desk (6.5"-10")		[SPARM014]	£189
Tubular Flip Up ³		[SPARM015]	N/C
Omit Armrests		[SPARM007]	N/C

¹ Not available with Frame Depths less than 14"

² Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.

³ Not available with Arm Pads or Half-Folding Back

Arm Pads

Standard Arm Pad [SPARM026]	N/C
Waterfall Arm Pad [SPARM020]	£25
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i>	
Foam Grip [SPARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guards

Not available with Height Adjustable T-Arm Armrests.

Adult Composite Side Guards (14" - 16") [SPSDG002]	N/C
Pediatric Composite Side Guard (10" - 14") [SPSDG001]	N/C

Pelvic Positioning Belts


1" Auto Buckle Pelvic Positioning Belt [SPPOB004]	£68
1" Padded Auto Buckle Pelvic Positioning Belt [SPPOB005]	£151
1 1/2" Auto Buckle Pelvic Positioning Belt [SPPOB001]	£68
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [SPPOB002]	£151
2" Airline Buckle Pelvic Positioning Belt [SPPOB003]	£110
Small Flex-Mount Positioning Belt (16-24") [SPPOB007]	£188
Medium Flex-Mount Positioning Belt (19-29") [SPPOB008]	£188
Large Flex-Mount Positioning Belt (25-37") [SPPOB009]	£188
Belt Mounting Band Clamp [SPACC054]	£103

Anti-tips

Rear Anti-Tips [SPATP001]	£212
Omit [SPATP000]	N/C

Other Accessories

Tool Kit [SPACC005]	£20
Removable Underseat Pouch [SPACC015]	£49
<i>Not available with 10" - 13" Seat Widths.</i>	
Pediatric Ki Mobility Backpack [SPACC026]	£100
Ki Mobility Backpack [SPACC009]	£100
Black Neoprene Impact Guards [SPIMG001]	£86
Cane and Crutch Holder [SPACC007]	£91
Canopy with Hardware [SPACC042]	£272
<i>See image to the left</i>	



Frame Color

Options	Price	Colors																								
		Gloss														Matte						Meta				
		Black	Yellow	Light Pink	Dark Pink	Glacier Blue	Electric Green	Radiant Orange	Opalescent Violet	Purple Plum	Rich Teal	Metallic Orange	Granny Smith	Candy Red	Candy Blue	Candy Purple	Shadow Black	Tactical Green	Dry Earth Tan	White	Shadow Black	Bright Red	Royal Blue	Charcoal Gray	Iridescent Navy	Iridescent Burgundy
Frame	N/C																									
Mag Color																										
Color Accent ¹ Package	N/C																									

¹ Includes axle plates and forks. Universal axle plate only.

Hydrographics Base Color

Not available with Frame Color. Hydrographics Base Color required with Hydrographics. Zebra not available with Black.

Options	Price	Hydrographics					Hydrographics Base Color						
		AOR1 (Desert Digital Camo)	NEXT Camo	Pink Camo	Carbon Fiber	Muddy Girl	Zebra	Black	Blue	Granny Smith	Pink	Purple	Red
Frame	£725												
Base Color	N/C												



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Spark Front Frame Height Matrix																							
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork
14.5	13.5	Std	4"	15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"
15	14			15.5	14			16.5	15.5			3/4	17.5			16.5	3/4			18.5	17		
15.5	14	3/4	4"	15.5 ¹	14 ¹	3/4	4"	17.5	16	1 1/2	5"	18.5	17	1 1/2	6"	18.5	17.5	3/4	6"	20.5	n/a	1 1/2	7"
16	14.5			16	14.5			16	14.5			Std	17			15.5	Std			17.5	16		
16	15	1 1/2	4"	16.5 ¹	15 ¹	1 1/2	5"	16.5	15	3/4	5"		17.5	16	Std	6"		19.5	18	1 1/2	7"	19.5	18
16.5	15.5			17	15.5			16.5	15.5			3/4	17.5	16			1 1/2	6"	18			16.5	3/4
15.5	14.5	Std	5"	15.5	14	Std	5"	17	16	3/4	5"		17.5	16.5	1 1/2	6"			18.5	17.5	Std	7"	
16	15			16	14.5			17.5	16			1 1/2	18	16.5			Std	18	17	3/4			6"
16.5	15	3/4	5"	16	15	3/4	5"	18	16.5	1 1/2	6"		18.5	17.5	Std	7"		19.5	18		Std	7"	
17	15.5			16.5	15.5			16.5	15.5			1 1/2	17	15.5			Std	18.5	17	1 1/2			7"
17	16	1 1/2	6"	17	15.5	1 1/2	6"	17	15.5	Std	19		17.5	1 1/2	7"	19		18	3/4		7"	21	
17.5	16.5			17.5	16			17.5	16		3/4	18	16.5			3/4	6"	18		16.5		Std	7"
16.5	15.5	Std	6"	16.5	15	Std	6"	17.5	16	3/4		6"	18	16.5	Std			7"	20	19	3/4		
17	16			17	15.5			17	15.5		3/4		18	16.5		3/4	6"		18.5	17		Std	7"
17.5	16	3/4	6"	17	16	3/4	6"	18.5	17	1 1/2		6"	19	17.5	Std			7"	20	18.5	1 1/2		
18	16.5			17.5	16.5			17.5	16.5		1 1/2		18	16.5		1 1/2	6"		18	16.5		1 1/2	7"
18	17	1 1/2	6"	18	16.5	1 1/2	6"	18.5	17	1 1/2		6"	18.5	17.5	1 1/2			7"	20.5	19	1 1/2		
18.5	17.5			18.5	17			18.5	17		1 1/2		19	17.5		1 1/2	6"		19	18		3/4	7"
		18	16.5	18	16.5	Std	18	16.5	Std	19.5		18.5	1 1/2	7"						7"			
		18.5	17	18.5	17		3/4	18.5		17	Std	20			19	1 1/2	7"						7"
		18.5	17.5	19	17.5	3/4		19	17.5	3/4		7"	19.5	18	1 1/2			7"					
		19	18	18.5	17.5		1 1/2	19	18		3/4		7"	20		18.5	1 1/2		7"				
		19.5	18.5	19	18	1 1/2		19	18	3/4		7"		20.5	19	1 1/2		7"					
		20	19	19.5	18.5		1 1/2	19.5	18		1 1/2		7"	21	n/a		1 1/2		7"				
				19.5	18	1 1/2		20	18.5	1 1/2		7"						7"					
				20	18.5		1 1/2	20.5	19		1 1/2		7"						7"				
				20.5	19	1 1/2						7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"						7"				
												7"						7"					
													7"										



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Axiom for Kids P - Positioning Cushion

Weight Capacity 165 lbs / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	£300									
Depth	8"	PXP0808BLK	PXP0908BLK	PXP1008BLK						
		PXP0808BLU	PXP0908BLU	PXP1008BLU						
		PXP0808PUR	PXP0908PUR	PXP1008PUR						
	9"	PXP0809BLK	PXP0909BLK	PXP1009BLK	PXP1109BLK					
		PXP0809BLU	PXP0909BLU	PXP1009BLU	PXP1109BLU					
		PXP0809PUR	PXP0909PUR	PXP1009PUR	PXP1109PUR					
	10"	PXP0810BLK	PXP0910BLK	PXP1010BLK	PXP1110BLK	PXP1210BLK	PXP1310BLK			
		PXP0810BLU	PXP0910BLU	PXP1010BLU	PXP1110BLU	PXP1210BLU	PXP1310BLU			
		PXP0810PUR	PXP0910PUR	PXP1010PUR	PXP1110PUR	PXP1210PUR	PXP1310PUR			
	11"			PXP1011BLK	PXP1111BLK	PXP1211BLK	PXP1311BLK	PXP1411BLK		
				PXP1011BLU	PXP1111BLU	PXP1211BLU	PXP1311BLU	PXP1411BLU		
				PXP1011PUR	PXP1111PUR	PXP1211PUR	PXP1311PUR	PXP1411PUR		
	12"			PXP1012BLK	PXP1112BLK	PXP1212BLK	PXP1312BLK	PXP1412BLK		
				PXP1012BLU	PXP1112BLU	PXP1212BLU	PXP1312BLU	PXP1412BLU		
				PXP1012PUR	PXP1112PUR	PXP1212PUR	PXP1312PUR	PXP1412PUR		
	13"			PXP1013BLK	PXP1113BLK	PXP1213BLK	PXP1313BLK	PXP1413BLK	PXP1513BLK	
				PXP1013BLU	PXP1113BLU	PXP1213BLU	PXP1313BLU	PXP1413BLU	PXP1513BLU	
				PXP1013PUR	PXP1113PUR	PXP1213PUR	PXP1313PUR	PXP1413PUR	PXP1513PUR	
	14"			PXP1014BLK	PXP1114BLK	PXP1214BLK	PXP1314BLK	PXP1414BLK	PXP1514BLK	PXP1614BLK
				PXP1014BLU	PXP1114BLU	PXP1214BLU	PXP1314BLU	PXP1414BLU	PXP1514BLU	PXP1614BLU
				PXP1014PUR	PXP1114PUR	PXP1214PUR	PXP1314PUR	PXP1414PUR	PXP1514PUR	PXP1614PUR
	15"					PXP1215BLK	PXP1315BLK	PXP1415BLK	PXP1515BLK	PXP1615BLK
						PXP1215BLU	PXP1315BLU	PXP1415BLU	PXP1515BLU	PXP1615BLU
						PXP1215PUR	PXP1315PUR	PXP1415PUR	PXP1515PUR	PXP1615PUR
16"					PXP1216BLK	PXP1316BLK	PXP1416BLK	PXP1516BLK	PXP1616BLK	
					PXP1216BLU	PXP1316BLU	PXP1416BLU	PXP1516BLU	PXP1616BLU	
					PXP1216PUR	PXP1316PUR	PXP1416PUR	PXP1516PUR	PXP1616PUR	
17"						PXP1317BLK	PXP1417BLK	PXP1517BLK	PXP1617BLK	
						PXP1317BLU	PXP1417BLU	PXP1517BLU	PXP1617BLU	
						PXP1317PUR	PXP1417PUR	PXP1517PUR	PXP1617PUR	
18"						PXP1318BLK	PXP1418BLK	PXP1518BLK	PXP1618BLK	
						PXP1318BLU	PXP1418BLU	PXP1518BLU	PXP1618BLU	
						PXP1318PUR	PXP1418PUR	PXP1518PUR	PXP1618PUR	



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Axiom for Kids SP - Skin Protection and Positioning Cushion

Weight Capacity 165 lbs. / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	£425									
Depth	8"	PXSP0808BLK	PXSP0908BLK	PXSP1008BLK						
		PXSP0808BLU	PXSP0908BLU	PXSP1008BLU						
		PXSP0808PUR	PXSP0908PUR	PXSP1008PUR						
	9"	PXSP0809BLK	PXSP0909BLK	PXSP1009BLK	PXSP1109BLK					
		PXSP0809BLU	PXSP0909BLU	PXSP1009BLU	PXSP1109BLU					
		PXSP0809PUR	PXSP0909PUR	PXSP1009PUR	PXSP1109PUR					
	10"	PXSP0810BLK	PXSP0910BLK	PXSP1010BLK	PXSP1110BLK	PXSP1210BLK	PXSP1310BLK			
		PXSP0810BLU	PXSP0910BLU	PXSP1010BLU	PXSP1110BLU	PXSP1210BLU	PXSP1310BLU			
		PXSP0810PUR	PXSP0910PUR	PXSP1010PUR	PXSP1110PUR	PXSP1210PUR	PXSP1310PUR			
	11"			PXSP1011BLK	PXSP1111BLK	PXSP1211BLK	PXSP1311BLK	PXSP1411BLK		
				PXSP1011BLU	PXSP1111BLU	PXSP1211BLU	PXSP1311BLU	PXSP1411BLU		
				PXSP1011PUR	PXSP1111PUR	PXSP1211PUR	PXSP1311PUR	PXSP1411PUR		
	12"			PXSP1012BLK	PXSP1112BLK	PXSP1212BLK	PXSP1312BLK	PXSP1412BLK		
				PXSP1012BLU	PXSP1112BLU	PXSP1212BLU	PXSP1312BLU	PXSP1412BLU		
				PXSP1012PUR	PXSP1112PUR	PXSP1212PUR	PXSP1312PUR	PXSP1412PUR		
	13"			PXSP1013BLK	PXSP1113BLK	PXSP1213BLK	PXSP1313BLK	PXSP1413BLK	PXSP1513BLK	
				PXSP1013BLU	PXSP1113BLU	PXSP1213BLU	PXSP1313BLU	PXSP1413BLU	PXSP1513BLU	
				PXSP1013PUR	PXSP1113PUR	PXSP1213PUR	PXSP1313PUR	PXSP1413PUR	PXSP1513PUR	
	14"			PXSP1014BLK	PXSP1114BLK	PXSP1214BLK	PXSP1314BLK	PXSP1414BLK	PXSP1514BLK	PXSP1614BLK
				PXSP1014BLU	PXSP1114BLU	PXSP1214BLU	PXSP1314BLU	PXSP1414BLU	PXSP1514BLU	PXSP1614BLU
				PXSP1014PUR	PXSP1114PUR	PXSP1214PUR	PXSP1314PUR	PXSP1414PUR	PXSP1514PUR	PXSP1614PUR
	15"					PXSP1215BLK	PXSP1315BLK	PXSP1415BLK	PXSP1515BLK	PXSP1615BLK
						PXSP1215BLU	PXSP1315BLU	PXSP1415BLU	PXSP1515BLU	PXSP1615BLU
						PXSP1215PUR	PXSP1315PUR	PXSP1415PUR	PXSP1515PUR	PXSP1615PUR
	16"					PXSP1216BLK	PXSP1316BLK	PXSP1416BLK	PXSP1516BLK	PXSP1616BLK
						PXSP1216BLU	PXSP1316BLU	PXSP1416BLU	PXSP1516BLU	PXSP1616BLU
						PXSP1216PUR	PXSP1316PUR	PXSP1416PUR	PXSP1516PUR	PXSP1616PUR
	17"						PXSP1317BLK	PXSP1417BLK	PXSP1517BLK	PXSP1617BLK
							PXSP1317BLU	PXSP1417BLU	PXSP1517BLU	PXSP1617BLU
							PXSP1317PUR	PXSP1417PUR	PXSP1517PUR	PXSP1617PUR
	18"						PXSP1318BLK	PXSP1418BLK	PXSP1518BLK	PXSP1618BLK
							PXSP1318BLU	PXSP1418BLU	PXSP1518BLU	PXSP1618BLU
							PXSP1318PUR	PXSP1418PUR	PXSP1518PUR	PXSP1618PUR

Accessories

Extra Cover

Extra Cover Black	£85
Extra Cover Blue	£85
Extra Cover Purple	£85

Solid Seat Insert


Solid Seat Insert	£85
-------------------	-----

Axiom for Kids Back Width (all models)	Fits back Frame Widths
10" (25.5 cm)	9" - 11" (23 cm - 28 cm)
12" (30.5 cm)	11" - 13" (28 cm - 33 cm)
14" (35.5 cm)	13" - 15" (33 cm - 38 cm)
16" (40.5 cm)	15" - 17" (38 cm - 43 cm)



Axiom for Kids PL - Posterior Lateral Back

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Trunk Widths (F):	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Contour Depth (H):	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)




Back

Includes Lumbar Pad. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Width			
		10	12	14	16
Retail		£550			
Height	10	PXPL1010			
	12	PXPL1012	PXPL1212	PXPL1412	
	14	PXPL1014	PXPL1214	PXPL1414	PXPL1614
	16		PXPL1216	PXPL1416	PXPL1616
	18				PXPL1618

Axiom for Kids AL - Posterior Back for Adjustable Laterals

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	8" (20cm)	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)
Component Lateral Width Range (J)*: <i>(*-1" with included spacers)</i>	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contour Depth (H):	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)



Back

Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Width			
		10"	12"	14"	16"
Retail		£550			
Height	10"	PXAL1010			
	12"	PXAL1012	PXAL1212	PXAL1412	
	14"	PXAL1014	PXAL1214	PXAL1414	PXAL1614
	16"		PXAL1216	PXAL1416	PXAL1616
	18"				PXAL1618

Axiom for Kids - Adjustable Laterals Trunk Supports (for PXAL Back)

Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)

Description	Item Code	Retail
3 x 4 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair)	PLF3x4PP	£250
4 x 5 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair)	PLF4x5PP	£250
3 x 4 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair)	PLTAS3x4PP	£500
4 x 5 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair)	PLTAS4x5PP	£500
3 x 4 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair)	PLTASC3x4PP	£500
4 x 5 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair)	PLTASC4x5PP	£500


Axiom for Kids - Adjustable Lateral Pelvic Supports (for PXAL Back)

Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)

Description	Item Code	Retail
3 x 4 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair)	PLPF3x4PP	£250
4 x 5 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair)	PLPF4x5PP	£250
3 x 4 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair)	PLASP3x4PP	£500
4 x 5 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair)	PLASP4x5PP	£500

Axiom for Kids ALCL - Posterior Back for Contoured Adjustable Laterals

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contoured Lateral Width Range (J):	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contour Depth (K):	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)



Back

Requires Contoured Adjustable Lateral Trunk Support selection below. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Width			
		10"	12"	14"	16"
Retail		£550			
Height	10"	PXAL1010CL			
	12"	PXAL1012CL	PXAL1212CL	PXAL1412CL	
	14"	PXAL1014CL	PXAL1214CL	PXAL1414CL	PXAL1614CL
	16"		PXAL1216CL	PXAL1416CL	PXAL1616CL
	18"				PXAL1618CL

Axiom for Kids - Contoured Adjustable Lateral Trunk Supports (for PXALCL Back)

Description	Item Code	Retail
3" Contoured Lateral Wings with Fixed Mounting Hardware (Pair)	PCL3	£250
4" Contoured Lateral Wings with Fixed Mounting Hardware (Pair)	PCL4	£250

4" Contoured Lateral Wings are not available on 10" or 12" Back Heights

Accessories

Axiom for Kids Back Accessories & Modifications	
Lumbar Pad	£40
Privacy Flap	£40
<i>Privacy Flap not available for 10" wide backs</i>	
Axiom DEX for Kids Hardware Extension	N/C
<i>Hardware extension adds an additional 2" of anterior depth adjustment</i>	
Axiom for Kids Adjustable Shoulder Guides	£50

Headrests								
Description	Standard Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		Posalinc Flip-Back Adjustable Hardware	
	Item Code	Retail	Item Code	Retail	Item Code	Retail	Item Code	Retail
Comfort Soft Headrest 8"	CS08SS	£268	CS08MA	£268	CS08PL	£268	CS08PLF	£268
Comfort Soft Headrest 10"	CS10SS	£268	CS10MA	£268	CS10PL	£268	CS10PLF	£268
Comfort Soft Headrest 14"	CS14SS	£268	CS14MA	£268	CS14PL	£268	CS14PLF	£268
Contour Positioning Pad, 8" X 6.25"	CP08SS	£303	CP08MA	£303	CP08PL	£303	CP08PLF	£303
Contour Positioning Pad, 9" X 6.5"	CP09SS	£303	CP09MA	£303	CP09PL	£303	CP09PLF	£303
Contour Positioning Pad, 10" X 7.5"	CP10SS	£303	CP10MA	£303	CP10PL	£303	CP10PLF	£303
Tri-soft Headrest 10"	TS10SS	£386	TS10MA	£386	TS10PL	£386	TS10PLF	£386
Tri-soft Headrest 14"	TS14SS	£386	TS14MA	£386	TS14PL	£386	TS14PLF	£386

Headrest Hardware		
<i>Not required with Axiom Headrests</i>		
Description	Item Code	Retail
Universal Headrest Mounting Plate	UMP	£43

Anterior Trunk Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	14941	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	14942	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	14943	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	14924	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	14927	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	14928	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	14851	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	14852	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	14853	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	14854	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	14856	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	14857	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	14858	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	14859	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	14282	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	14283	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	14284	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	14257	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	14258	£300



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com



Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Opalescent Violet Back Shell	N/C
Radiant Orange Back Shell	N/C
Electric Green Back Shell	N/C
Glacier Blue Back Shell	N/C
Dark Pink Back Shell	N/C
Light Pink Back Shell	N/C
Rich Teal Back Shell	N/C
Purple Plum Back Shell	N/C
Black Back Shell	N/C
Candy Blue Back Shell	N/C
Candy Purple Back Shell	N/C
Candy Red Back Shell	N/C
Granny Smith Back Shell	N/C
Metallic Orange Back Shell	N/C
Yellow Back Shell	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Bright Red Back Shell	N/C
White Back Shell	N/C
Dry Earth Tan Back Shell	N/C
Tactical Green Back Shell	N/C
Charcoal Gray Back Shell	N/C
Iridescent Burgundy Back Shell	N/C
Iridescent Navy Back Shell	N/C
Royal Blue Back Shell	N/C

Metallic Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Blue Steel Back Shell	N/C
-----------------------	-----