



SPARKY

The Sparky® is a rear-wheel drive power chair designed to deliver rehab capability in a compact, powerful paediatric package. Engineered to operate in tight spaces and to keep up with today's active child. Standard features include powered seat tilt and attendant operated manual recline.

QUANTUM®

Specifications

WEIGHT CAPACITY ⁴	8 stone 13 lbs. (56.7 kg)
MAXIMUM SPEED ¹	Up to 4 mph
GROUND CLEARANCE ²	95 mm (3.75")
TURNING RADIUS ²	743 mm (29.25")
OVERALL SIZE ^{2, 6}	Length: 959 mm (37.75") with front riggings Width: 610 mm (24")
DRIVE WHEELS	304 mm (12") pneumatic, standard
CASTER WHEELS	Front: 203 mm (8") pneumatic, standard
ANTI-TIP WHEELS	76 mm (3")
SUSPENSION	Rear-wheel drive
COMPONENT WEIGHTS ^{5, 7}	Base: 130.5 lbs (59.2 kg) Batteries: 23 lbs (10.5 kg) each
DRIVETRAIN	Two motor, Rear-Wheel
BRAKES	"Intelligent Braking" electronic regenerative, disc park brake
BATTERIES ³	Two 40 Ah
BATTERY CHARGER	5A, Off-board
MAXIMUM OBSTACLE CLIMBING ABILITY	40 mm (1.57")
SEAT	Dimensions: 330 mm (13") height of back canes 254-355 mm (10"-14") adjustable depth 254-355 mm (10"-14") adjustable width 432-476 mm (17"-18.75") seat to floor height

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability

Features

- ION Power Tilt
- Attendant Recline
- Crash tested
- Swing-away Joystick Mount
- Tension Adjustable Backrest
- Headrest
- VR2 Controls

Options

- ION Seat Sizes: 10" x 10", 12" x 12" and 14" x 14" (all width and depth adjustable)
- ADI Contoured Backrests
- Thoracic Supports
- iFit comfort harness
- iFit pelvic support
- Dual Controls
- Manual Elevating Legrests

Colours



White



Blue

* Please contact your representative for the full range of options and colours.