

32 Wedgwood Road,
Bicester, OX26 4UL

01869 324600

eu stealthproducts.com



EXCLUSIVE TO STEALTH PRODUCTS

No matter the season or the situation, the patented award-winning CoolCore fabric wicks moisture away from your body, helping keep you dry and comfortable. The chemical-free blends of the fabric also help with moisture transportation and moisture evaporation, making our cushions cool better and last longer.

STEALTH PRODUCTS CUSHION SPECIFICATIONS

Cushion Name	Weight Capacity	Weight of Cushion	Cushion Height	Width Range	Depth Range	RRP
Simplicity G	400 lbs	2 lbs	2"	10" - 24"	10" - 24"	£70 - £80
Zen SP	450 lbs	3.5 lbs	3"	10" - 24"	10" - 24"	£140 - £160
Glacial SP	400 lbs	2.5 lbs	2.5"	10" - 24"	10" - 24"	£180 - £200
Premiere P	550 lbs	3.5 lbs	3.5"	16" - 22"	16" - 22"	£200 - £220
Spectrum Foam SPP	650 lbs	3 lbs	3"	16" - 32"	16" - 28"	£250 - £300
Spectrum Gel SPP	550 lbs	3.5 lbs	3"	14" - 24"	16" - 24"	£250 - £300
Solution SPP	500 lbs	5.4 lbs	3"	10" - 24"	10" - 22"	£180 - £200
Tru Comfort 2 SPP	450 lbs	4 lbs	4.5"	16" - 22"	16" - 24"	£150 - £170
Essence SPP	500 lbs	3.5 - 5 lbs	3.5"	16" - 22"	16" - 24"	£200 - £230

**SOLUTIONS
FOR LIFE**





CUSHION COLOUR GUIDE

To make your cushion selection process easier, Stealth Products Cushions are marked with a colour coded patch. The colour of the patch identifies the function of the cushion.

Cushion RRP:

- Simplicity G:** £70 - £80
- Zen SP:** £140 - £160
- Glacial SP:** £180 - £200
- Premier P:** £200 - £220
- Spectrum Foam SPP:** £250 - £300
- Spectrum Gel SPP:** £250 - £300
- Solution SPP:** £180 - £200
- Tru-Comfort 2 SPP:** £150 - £170
- Essence SPP:** £200 - £230



eu.stealthproducts.com



SIMPLICITY G

Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 10" - 24"W x 10" - 24"D



ZEN SP

Features:

- Slightly Contoured, high-density molded foam base
- FRG Comfort foam topper for skin protection
- Standard Visco topper (exclusive FRG SunMate) for added comfort
- FRG SunMate is an all natural flame retardant material.
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- Optional Urethane liner
- 4-Way Stretch Poly Liner optional
- Sizes available - 10"-24"W x 10"-24"D
- Optional Mesh and Stealth Tek cover
- Optional ABS or ACM solid seat insert
- Special Order configurations available



GLACIAL SP

Features:

- Contoured, high-density molded foam base
- Soft gel top layer for skin protection
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- 4-Way Stretch Poly or Urethane liner optional
- Sizes available - 10"-24" W x 10"-24"D
- Optional Mesh and Stealth Tek cover
- Optional ABS or ACM solid seat insert



PREMIERE P

Features:

- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base
- Built-in positioning platform for off loading
- Standard Visco topper (exclusive FRG SunMate) for added comfort
- FRG SunMate is an all natural flame retardant material.
- Coccyx and sacral relief
- Sizes available - 16"-22"W x 16"-22"D
- CoolCore® technology cover standard.
- Optional Urethane liner
- Optional 4-Way Stretch Poly liner
- Optional Mesh outer cover
- Optional ABS or ACM Solid Seat Insert



SPECTRUM FOAM SPP

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 16" - 32"W x 16" - 28"D
- Custom sizes available
- Optional Stealth-Tek, Mesh cover



SPECTRUM GEL SPP

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 14" - 24"W x 16" - 24"D
- Custom sizes available
- Optional Stealth-Tek, Mesh cover



SOLUTION SPP

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 22"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.
- 4-Way Stretch Poly liner optional



TRU-COMFORT 2 SPP

Features:

- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 16" - 24"D
- Custom sizes available
- 4-Way Stretch Poly liner optional



ESSENCE SPP

Features:

- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base.
- Coccyx and sacral relief.
- CoolCore® technology cover standard.
- 4-Way Stretch Poly liner optional.
- Sizes available - 16"-22"W x 16"-24"D
- Special Order configurations available
- Adjustable contour wedge package provides better positioning customization for the pelvic region. Wedge height from 0.5" - 2" with the wedge kit.
- Optional mesh cover.
- Solid Seat insert included.