

Start M6 Junior

The Start M6 Junior is a folding basic wheelchair for children who only depend on a wheelchair for a limited amount of time. It is robust and equally well suited for short distances outdoors and in the home.

Start M6 Junior Standard Model

- Upholstery colour black
- Aluminium tube frame
- Nylon seat and back upholstery, adjustable back
- Side panel with height adjustable arm rest
- Rear wheels with quick release axles
- Knee lever wheel lock
- Swing-away footrest holders
- Plastic footplate
- Individual footrest, angle adjustable
- Tip-assist integrated into the frame



Start M6 Junior

Options

- Headrest
- Side thoracic support (right/left)
- Back angle adjustment 30°
- Push handles, various models
- Tray
- Wheel lock lever extension
- Push rings, various models
- Spoke protectors
- Footrest, various models
- Anti-tipper (plug-on, swing-away)



Attendant brake option



Push handle options



Short lower leg length footplate



Swing away footrest

Technical Data

	Start M6 Junior
Max user weight	90kg
Seat width	280mm – 380mm
Seat depth	310mm – 370mm
Back height	300mm – 375mm
Seat height	Front 400 – 490mm Rear 370mm – 450mm
Lower leg length	160mm – 450mm
Variable armrest height	Up to 400mm