

















# PDG MOBILITY PRODUCTS

<i>ELEVATION</i>	<i>FUZE T20</i>	<i>FUZE T50</i>	<i>ECLIPSE</i>	<i>STELLAR</i>	<i>STELLAR GL</i>	<i>STELLAR LEAP</i>	<i>BENTLEY</i>																																																																
																																																																							
																																																																							
<table border="0"> <tr> <td>SEAT WIDTH 12"-19"</td> <td>DYNAMIC SEAT HEIGHT Up to 10" of range</td> </tr> <tr> <td>MAXIMUM CAPACITY 250 lbs/113 kg</td> <td>FRONT SEAT HEIGHT 18"-21"</td> </tr> <tr> <td>BACK ANGLE 85°-115°</td> <td>WHEEL CAMBER 0°, 2°, 4°, 6°</td> </tr> <tr> <td>CHAIR WEIGHT INC WHEELS 22 lbs/10 kg</td> <td>SEAT DEPTH 14"-20"</td> </tr> </table>	SEAT WIDTH 12"-19"	DYNAMIC SEAT HEIGHT Up to 10" of range	MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 18"-21"	BACK ANGLE 85°-115°	WHEEL CAMBER 0°, 2°, 4°, 6°	CHAIR WEIGHT INC WHEELS 22 lbs/10 kg	SEAT DEPTH 14"-20"	<table border="0"> <tr> <td>SEAT WIDTH 12"-22"</td> <td>DYNAMIC TILT Up to 20°</td> </tr> <tr> <td>MAXIMUM CAPACITY 250 lbs/113 kg</td> <td>FRONT SEAT HEIGHT 13"-20"</td> </tr> <tr> <td>BACK ANGLE 84°-120°</td> <td>OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")</td> </tr> <tr> <td>SEAT DEPTH 15"-20"</td> <td>CHAIR WEIGHT 42 lbs/19 kg</td> </tr> </table>	SEAT WIDTH 12"-22"	DYNAMIC TILT Up to 20°	MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 13"-20"	BACK ANGLE 84°-120°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")	SEAT DEPTH 15"-20"	CHAIR WEIGHT 42 lbs/19 kg	<table border="0"> <tr> <td>ADJ SEAT WIDTH 12"-22" FIXED SEAT WIDTH 15"-24"</td> <td>DYNAMIC TILT Up to 50°</td> </tr> <tr> <td>MAX CAPACITY ADJ FRAME 250 lbs/113 kg MAX CAPACITY FIXED FRAME 350 lbs/158 kg</td> <td>FRONT SEAT HEIGHT 13"-20"</td> </tr> <tr> <td>BACK ANGLE 84°-120°</td> <td>OVERALL WIDTH +7.5" (RW 12/16) +8.5" (RW 20")</td> </tr> <tr> <td>SEAT DEPTH 15"-22" ALSO AVAILABLE 23" &amp; 24"</td> <td>CHAIR WEIGHT 46 lbs/21 kg</td> </tr> </table>	ADJ SEAT WIDTH 12"-22" FIXED SEAT WIDTH 15"-24"	DYNAMIC TILT Up to 50°	MAX CAPACITY ADJ FRAME 250 lbs/113 kg MAX CAPACITY FIXED FRAME 350 lbs/158 kg	FRONT SEAT HEIGHT 13"-20"	BACK ANGLE 84°-120°	OVERALL WIDTH +7.5" (RW 12/16) +8.5" (RW 20")	SEAT DEPTH 15"-22" ALSO AVAILABLE 23" & 24"	CHAIR WEIGHT 46 lbs/21 kg	<table border="0"> <tr> <td>SEAT WIDTH 20"-48"</td> <td>BACK HEIGHT 25" &amp; 30"</td> </tr> <tr> <td>MAX CAPACITY ADJ FRAME 600 lbs/272 kg HD PACKAGE 1000 lbs/453 kg</td> <td>FRONT SEAT HEIGHT 14.5"-20"</td> </tr> <tr> <td>BACK ANGLE 85°-100°</td> <td>OVERALL WIDTH Seat Width +7.5"</td> </tr> <tr> <td>SEAT DEPTH 16"-24"</td> <td>CHAIR WEIGHT 44 lbs/20 kg</td> </tr> </table>	SEAT WIDTH 20"-48"	BACK HEIGHT 25" & 30"	MAX CAPACITY ADJ FRAME 600 lbs/272 kg HD PACKAGE 1000 lbs/453 kg	FRONT SEAT HEIGHT 14.5"-20"	BACK ANGLE 85°-100°	OVERALL WIDTH Seat Width +7.5"	SEAT DEPTH 16"-24"	CHAIR WEIGHT 44 lbs/20 kg	<table border="0"> <tr> <td>SEAT WIDTH 14"-32"</td> <td>DYNAMIC TILT Up to 45°</td> </tr> <tr> <td>MAX CAPACITY ADJ FRAME 250 lbs/113 kg STELLAR HD 600 lbs/272 kg</td> <td>FRONT SEAT HEIGHT 14"-20"</td> </tr> <tr> <td>BACK ANGLE 90°-122°</td> <td>OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20") STELLAR HD +11.1" (RW 24)</td> </tr> <tr> <td>SEAT DEPTH 16"-22"</td> <td>CHAIR WEIGHT 53 lbs/24 kg</td> </tr> </table>	SEAT WIDTH 14"-32"	DYNAMIC TILT Up to 45°	MAX CAPACITY ADJ FRAME 250 lbs/113 kg STELLAR HD 600 lbs/272 kg	FRONT SEAT HEIGHT 14"-20"	BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20") STELLAR HD +11.1" (RW 24)	SEAT DEPTH 16"-22"	CHAIR WEIGHT 53 lbs/24 kg	<table border="0"> <tr> <td>SEAT WIDTH 14"-22"</td> <td>DYNAMIC TILT Up to 45°</td> </tr> <tr> <td>MAXIMUM CAPACITY 250 lbs/113 kg</td> <td>FRONT SEAT HEIGHT 13"-18"</td> </tr> <tr> <td>BACK ANGLE 90°-122°</td> <td>OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")</td> </tr> <tr> <td>SEAT DEPTH 16"-22"</td> <td>CHAIR WEIGHT 50 lbs/23 kg</td> </tr> </table>	SEAT WIDTH 14"-22"	DYNAMIC TILT Up to 45°	MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 13"-18"	BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")	SEAT DEPTH 16"-22"	CHAIR WEIGHT 50 lbs/23 kg	<table border="0"> <tr> <td>SEAT WIDTH 14"-22"</td> <td>ANTERIOR TILT Up to 30° POSTERIOR TILT Up to 20°</td> </tr> <tr> <td>MAXIMUM CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg</td> <td>FRONT SEAT HEIGHT 13"-18"</td> </tr> <tr> <td>BACK ANGLE 90°-122°</td> <td>OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")</td> </tr> <tr> <td>SEAT DEPTH 16"-22"</td> <td>CHAIR WEIGHT 58 lbs/26 kg</td> </tr> </table>	SEAT WIDTH 14"-22"	ANTERIOR TILT Up to 30° POSTERIOR TILT Up to 20°	MAXIMUM CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg	FRONT SEAT HEIGHT 13"-18"	BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")	SEAT DEPTH 16"-22"	CHAIR WEIGHT 58 lbs/26 kg	<table border="0"> <tr> <td>SEAT WIDTH 14"-32"</td> <td>DYNAMIC TILT Up to 20°</td> </tr> <tr> <td>MAX CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg</td> <td>FRONT SEAT HEIGHT 13"-20"</td> </tr> <tr> <td>BACK ANGLE 90°-122°</td> <td>OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")</td> </tr> <tr> <td>SEAT DEPTH 16"-22"</td> <td>CHAIR WEIGHT 51 lbs/23 kg</td> </tr> </table>	SEAT WIDTH 14"-32"	DYNAMIC TILT Up to 20°	MAX CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg	FRONT SEAT HEIGHT 13"-20"	BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")	SEAT DEPTH 16"-22"	CHAIR WEIGHT 51 lbs/23 kg
SEAT WIDTH 12"-19"	DYNAMIC SEAT HEIGHT Up to 10" of range																																																																						
MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 18"-21"																																																																						
BACK ANGLE 85°-115°	WHEEL CAMBER 0°, 2°, 4°, 6°																																																																						
CHAIR WEIGHT INC WHEELS 22 lbs/10 kg	SEAT DEPTH 14"-20"																																																																						
SEAT WIDTH 12"-22"	DYNAMIC TILT Up to 20°																																																																						
MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 13"-20"																																																																						
BACK ANGLE 84°-120°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")																																																																						
SEAT DEPTH 15"-20"	CHAIR WEIGHT 42 lbs/19 kg																																																																						
ADJ SEAT WIDTH 12"-22" FIXED SEAT WIDTH 15"-24"	DYNAMIC TILT Up to 50°																																																																						
MAX CAPACITY ADJ FRAME 250 lbs/113 kg MAX CAPACITY FIXED FRAME 350 lbs/158 kg	FRONT SEAT HEIGHT 13"-20"																																																																						
BACK ANGLE 84°-120°	OVERALL WIDTH +7.5" (RW 12/16) +8.5" (RW 20")																																																																						
SEAT DEPTH 15"-22" ALSO AVAILABLE 23" & 24"	CHAIR WEIGHT 46 lbs/21 kg																																																																						
SEAT WIDTH 20"-48"	BACK HEIGHT 25" & 30"																																																																						
MAX CAPACITY ADJ FRAME 600 lbs/272 kg HD PACKAGE 1000 lbs/453 kg	FRONT SEAT HEIGHT 14.5"-20"																																																																						
BACK ANGLE 85°-100°	OVERALL WIDTH Seat Width +7.5"																																																																						
SEAT DEPTH 16"-24"	CHAIR WEIGHT 44 lbs/20 kg																																																																						
SEAT WIDTH 14"-32"	DYNAMIC TILT Up to 45°																																																																						
MAX CAPACITY ADJ FRAME 250 lbs/113 kg STELLAR HD 600 lbs/272 kg	FRONT SEAT HEIGHT 14"-20"																																																																						
BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20") STELLAR HD +11.1" (RW 24)																																																																						
SEAT DEPTH 16"-22"	CHAIR WEIGHT 53 lbs/24 kg																																																																						
SEAT WIDTH 14"-22"	DYNAMIC TILT Up to 45°																																																																						
MAXIMUM CAPACITY 250 lbs/113 kg	FRONT SEAT HEIGHT 13"-18"																																																																						
BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")																																																																						
SEAT DEPTH 16"-22"	CHAIR WEIGHT 50 lbs/23 kg																																																																						
SEAT WIDTH 14"-22"	ANTERIOR TILT Up to 30° POSTERIOR TILT Up to 20°																																																																						
MAXIMUM CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg	FRONT SEAT HEIGHT 13"-18"																																																																						
BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")																																																																						
SEAT DEPTH 16"-22"	CHAIR WEIGHT 58 lbs/26 kg																																																																						
SEAT WIDTH 14"-32"	DYNAMIC TILT Up to 20°																																																																						
MAX CAPACITY 250 lbs/113 kg HD PACKAGE 450 lbs/204 kg	FRONT SEAT HEIGHT 13"-20"																																																																						
BACK ANGLE 90°-122°	OVERALL WIDTH +9.5" (RW 12/16) +10.5" (RW 20")																																																																						
SEAT DEPTH 16"-22"	CHAIR WEIGHT 51 lbs/23 kg																																																																						
Anti-shear side guard system	20° of tilt without knee rise	Adjustable axle depth and seat height	Forward wheelbase for maneuverability	Extra wide heavy duty version available (Stellar HD)	Anterior and posterior setups available	30° of linked dynamic back recline	20° of tilt without knee rise																																																																
Dynamic seat-height system	Foot and handle propel in tilt position	Constant center of gravity	Removable back and rear wheels for easy transport	Knees remain low in tilt	Knees remain low in tilt	30° of anterior tilt with automatic back adjustment	Foot and hand propel in tilt position																																																																
On-the-fly adjustable back-angle	Adjustable seat width and seat depth	Power Tilt available	Dynamic back recline option	Low seat-to-floor height	Power Tilt available	Anterior lift force balanced to user's weight	Low seat-to-floor height																																																																
Fold down back	Low seat-to-floor height	Pediatric version available (Fuze JR)	Narrow overall width	Ultra-durable version available (Stellar Impact)	Durable construction for heavy use	Power Tilt available	Durable construction for heavy use																																																																
Ultralight aluminum construction	Transit tie-down	Transit tie-down	HD Package available for 1000 lbs capacity	Transit tie-down *(14" - 22" only)	Transit tie-down	Transit tie-down	Transit tie-down *(14" - 22" only)																																																																