


Tarta Kid

rev. 06- 08.03.19

Technical *sheet*

English

Visit our website
www.tartadesign.it
and discover our
products.





Main features

Design

We are convinced that aesthetics and technical characteristics are equally important. The "Tarta Kid" design has been meticulously studied to be essential and at the same time lively.

Multipurpose

Tarta Kid has been designed to be used in different situations during the day. From chair to stroller, from home couch to restaurant table.



Recommended age

This product is suitable for 8 month-old to 3-4 year-old children.



Washable at 40°

Adjustable


Tarta kid has been thought to be regulated by parents. Very easy adjustments allow the immediate product management for everyone.

Comfort

The foam and the microfiber fabric give the pads a unique feeling of pleasantness and softness.




Height and depth




- | "Closed" cm 19
- | "Open" cm 26

Height: the values are measured by the level of the seat to the underarm measurement.





Inclination



Inclination degrees: -5° +30°

The backrest can be tilted by tightening or loosening the screw between the seat and the backrest.




Width



The width of the Kid can be easily changed without the use of any tools.. Thanks to the use of a special aluminum, the width can be changed with a simple manually-exerted pressure.

Approximate values:

- chest:** 13 to 23 cm
 - pelvis:** 17 to 25 cm
- 



Accessories



Standard headrest



Occipital headrest



Mandibular headrest



Articulated bracket



Pelvic belt



Chest belt



Chest support




Abductor



5 point harness




Lumbar wedge



Seat wedge



Flexibility block system



Footrest rise cushion



Seat rise cushion




Wooden high chair



Features:

- 1 adjustable footrest in height
- 2 height from the ground base: 55cm
- 3 inclined base of 8°


1



2



3






Height adjustable highchair



Features:

- ① height adjustable
seat height from 42 to 67 cm
- ② tilting base from 8° to 22°
- ③ equipped with wheels
- ④ foldable
- ⑤ removable table
- ⑥ adjustable footrest in height



Stroller



Features:

- ① aluminum ultra light frame
- ② foldable frame
- ③ tilting: 3 positions
- ④ facing parent mode
- ⑤ adjustable footrest in height






Table seat



Features:

- ① light to carry
- ② safe and fast application




Ground base



Features:

- ① easy to carry
- ② useful to position the baby correctly on the floor



Aeronautic aluminum mouldable

ideal to follow the requirements over time.

This material allows a simple and immediate management of Tarta Kid for anyone.

Green fabric

ecological microfiber soft to the touch that ensures maximum washability and resistance in time.


The fabric can be washed in the washing machine up to 40 °.



Pillow microspheres


extra-comfortable allowing multiple adaptation solutions.

With gentle pressure, you can create temporary positions to stabilize the child's arm.



Questions?
Info?
Help?

Visit our website:
www.tartadesign.it
or contact us at:
info@tartadesign.it



Tarta Design srl

Via Cottonificio, 56
33037 Pasian di Prato
Udine - Italy
T +39 0432 722 533
F +39 0432 542 920

www.tartadesign.it
info@tartadesign.it