



Date	UK (English)	
User	Company stamp	
End user weight (kg)		Folding frame wheelchair
Language user manual	Responsible	£ 2771,00
		Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm		
Seat width (LS)	<input type="checkbox"/> 36	<input type="checkbox"/> 39 <input type="checkbox"/> 42 <input type="checkbox"/> 45 <input type="checkbox"/> 48
Seat depth (PS)	<input type="checkbox"/> 37,5	<input type="checkbox"/> 40 <input type="checkbox"/> 42,5 <input type="checkbox"/> 45 <input type="checkbox"/> 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. - calculated with angle backrest 90° to the seat		
Footplate distance (DP) (cm) ±0.5cm - adjustment steps 1.5cm	
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		
Setting	☆ Standard <input type="checkbox"/> Prudential <input type="checkbox"/> Active	
camber	☆ 2° <input type="checkbox"/> 0° <input type="checkbox"/> 4°	
Backrest Height (TS)	Std.Value 33,5(min 24.5-max 47) adjustment steps 1.5cm	
Backrest angle (to the seat)	☆ 90° <input type="checkbox"/> -6° <input type="checkbox"/> 6°	
Front height (STD 49) - Rear height (STD 43)		


Front frame	
<input checked="" type="checkbox"/> Aluminium standard frame Angle 100° 	<input type="checkbox"/> Short frame - Ideal with Seat depth (PS)37,5 - 40 <input type="checkbox"/> Long frame - Ideal with Seat depth (PS)42,5 - 45 - 47,5
<input type="checkbox"/> Aluminium frame Angle 90° - (max 5" caster) 	<input type="checkbox"/> Short frame - Ideal with Seat depth (PS)37,5 - 40 <input type="checkbox"/> Long frame - Ideal with Seat depth (PS)42,5 - 45 - 47,5

Frame colour		
<input checked="" type="checkbox"/> Opaque black	<input type="checkbox"/> Tekno anthracite	<input type="checkbox"/> Matt grey
<input type="checkbox"/> Pearl white	<input type="checkbox"/> Matt Midnight blue	<input type="checkbox"/> Matt red
<input type="checkbox"/> Matt orange		

Backrest
<input checked="" type="checkbox"/> Adjustable velcro straps and padded transpiring upholstery

Handle tube	
<input checked="" type="checkbox"/> Tubes with standard handle	
<input type="checkbox"/> Folding handles	 02958 260,00
Not available with: Push handles for companion, Brake handles for companion (Drum Brake)	

Casters		
<input checked="" type="checkbox"/> 4" Soft		
<input type="checkbox"/> 3" Sport	<i>Lightest configuration</i>	
<input type="checkbox"/> 5" Soft		
<input type="checkbox"/> 6" Soft	(Available only with 100° angle)	
<input type="checkbox"/> 6" Pneumatic	(Available only with 100° angle)	
<input type="checkbox"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824 130,00
<input type="checkbox"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398 130,00
<input type="checkbox"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943 130,00
<input type="checkbox"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518 130,00
<input type="checkbox"/> 5" Pneumatic		10844 65,00

Support front fork	
<input checked="" type="checkbox"/> Standard support front fork	
<input type="checkbox"/> Front fork support V Design 2.0	 06850 152,00
Exclude the installation of QUICK RELEASE FRONT FORKS	

Date	Order	User	Address
------	-------	------	---------

Colour set: Includes rear wheel hub and rim, rear wheel plates and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black	22" 24" 25" 26"		
○ Bright silver	22" 24"	00833	130,00
○ Bright lime	22" 24"	00833	130,00
○ Bright blue	22" 24"	00833	130,00
○ Bright red	22" 24"	00833	130,00
○ Bright fuchsia	22" 24"	00833	130,00
○ Pearl white	22" 24"	00833	130,00
○ Tekno anthracite	22" 24"	00833	130,00
○ Matt Midnight blue	22" 24"	00833	130,00
○ Matt grey	22" 24"	00833	130,00
○ Matt red	22" 24"	00833	130,00

Rear wheels

☆ 24" Ultra-light Progeo - Standard!			
□ 22" Ultra-light Progeo - Standard!			
□ 25" Ultra-light Progeo - Standard!	(colour: black)		
□ 26" Progeo - Standard!	(colour: black)		
○ 24" ULTRA-LIGHT PROGEO KEY SPORT	<i>Lightest configuration</i>	04383	902,00
○ 24" Ultra-light "Spinergy LX"		01685	902,00
○ 24" Ultra-light "Spinergy LXL"	□ Black hub and white spokes □ Hub and spokes black	07846	902,00
○ 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.415,00
○ 25" Ultra-light "Spinergy LX"		02101	902,00
○ 25" Ultra-light "Spinergy LXL"	□ Black hub and white spokes	07847	902,00
○ 26" Ultra-light "Spinergy LX"		02102	902,00
○ 24" Off-road (aluminium handrim black)		01392	231,00
○ 22" with brake handles for companion (drum brake) (aluminium handrim)	Available only with: camber 0° and Setting Standard	03385	371,00
○ 24" with brake handles for companion (drum brake) (aluminium handrim)	Available only with: camber 0° and Setting Standard	03384	371,00
○ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°	○ Right ○ Left	02103	413,00

Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	<i>Lightest configuration</i>	
□ High pressure profile tyre	(24")		
○ Marathon Plus	(24"-25")	with anti-puncture protection	01669 65,00
○ Anti-puncture - Colour Black	(24")		11083 130,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"			
○ Anti-puncture - Colour Grey	(22")		04181 101,00


Handrim

Wheel distance ☆ 3 cm □ 4 cm

☆ Aluminium	(22"-24"-25"-26")	<i>Lightest configuration</i>	
○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170 130,00
○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684 152,00
○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016 283,00
○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017 324,00
○ Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834 359,00
○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023 411,00
Oval countour with a rubber insert			
○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025 271,00
Ergonomic black profile			
○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026 295,00
Ergonomic profile with rubber coating (ideal for quadriplegics)			
○ SURGE LT	(25")	Oval countour with a rubber insert	06427 371,00
○ Quad	(24")		00693 206,00
○ Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050 165,00


Wheel locks


☆ Bent lever push to lock			
□ Straight lever push to lock			
□ Straight lever pull to lock			
○ Sport brake		06428	130,00
○ Brake with extension lever push to lock		01949	65,00
○ Brake with extension lever pull to lock		03372	65,00






Date	Order	User	Address
------	-------	------	---------


Footplate - See reference tables N°3


- Automatic Closure footplate rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External  11526
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 195,00

- One piece footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External  11527
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 195,00

- One piece aluminium footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal 2/3 External External  11528

- Plastic footrest rake adjustable   06861 89,00
 Minimum footplate distance 35cm. 2/3 Internal 2/3 External
 Need a stretch calf band.
 NOT AVAILABLE WITH FRONT FORK SUPPORT V DESIGN 2.0

- One piece Carbon fibre footplate rake adjustable, swing away, flip-up.
 Minimum footplate distance 35cm. Internal External  03291 283,00


- Footplate COMPACT, automatic Closure, rake and depth adjustable  09957 260,00
 Minimum footplate distance 35cm. Internal External




Mudguard

- Mudguard**
 - ABS mudguard
 - ABS Straight sideguard
 - Carbon fibre mudguard
 - Carbon fibre Straight sideguard 03946 324,00
 08918 324,00

- Removable version**
 - Removable ABS mudguard
 - Removable ABS fibre straight sideguard
 - Removable carbon fibre mudguard
 - Removable Carbon fibre straight sideguard 08921 130,00
 08921 130,00
 08919 448,00
 08920 448,00

Armrests





- Tubular armrests - Steel support with safety locking system
 (Swing away, removable and height adjustable armrest)
 Not available with anti-tippers if the rear seat height is less than 40cm  00694 260,00



- Removable armrest with ABS mudguard**
 - Desk armrests PROGEO (Removable and height adjustable)  09956 231,00
 - Removable armrest with ABS straight sideguard  08881 231,00
 - Removable armrest with carbon-fiber straight sideguard  08882 467,00

Accessories

- Pad for front frame protectors 07058 54,00
- Backrest-seat connection cloth 10241 18,00
- Removable storage bag (with quick coupling) 09811 130,00

Date	Order	User	Address
------	-------	------	---------

Optional 1			
<u>Only one optional devices for side can be chosen</u>			
<input type="radio"/> Transit wheels	Exclude the installation of : Anti-tippers - Tipping aid - Crutches support		00836 152,00
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) <input type="radio"/> Right <input type="radio"/> Left		00023 130,00 00024 130,00
<input type="radio"/> Tipping aid	<input type="radio"/> Right <input type="radio"/> Left		00025 76,00 01046 76,00
<input type="radio"/> Crutches support	<input type="radio"/> Right <input type="radio"/> Left		00838 76,00 00837 76,00

Optional 2			
<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947 182,00
<input type="radio"/> Quick release front forks	To be matched only with: Standard support front fork		00030 152,00
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048 195,00
<input type="radio"/> Safety pelvic belt			03972 65,00
<input type="radio"/> Harness			04103 130,00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013 41,00
<input type="radio"/> Swing away lateral support Not available with Posture backrest PHYSIO	<input type="radio"/> Right <input type="radio"/> Left		02097 260,00 02098 260,00
<input type="radio"/> Spoke guards Black	(22"-24")		01043 152,00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835 152,00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019 182,00
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106 182,00








Additional Kit			
<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378 513,00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357 902,00

Date	Order	User	Address
------------	-------------	------------	---------------








PHYSIO - Posture backrest (Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- AIR** Carbon shell
not available for seat width 48

	BACKREST HEIGHT	28	34	40		1.078,00	
Shell's shape	Lateral retaining wing A=high - R=regular						
7,5 cm 	A	R	A	R	A	R	
							
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

- POSTURE** Aluminum shell

	BACKREST HEIGHT	34	40	52		1.043,00	
Shell's shape	Lateral retaining wing R=regular						
7,5 cm 	R	R	R	R	R	R	
							
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Configuration DEEP with side support PAD (height adjustable)


- Alluminum lateral with padded, breathable upholstery.



	Small (h 13cm)		Large (h 20cm)
		<input type="checkbox"/> Right 10272	
		30,00	<input type="checkbox"/> Right 10273
			41,00
		<input type="checkbox"/> Left 10272	
		30,00	<input type="checkbox"/> Left 10273
			41,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



<input type="radio"/> Large (18x12)		05647 188,00
<input type="radio"/> Small (14x9)		05643 188,00
<input type="radio"/> Lateral support extension	Quantity (number pcs).....	
<input type="radio"/> Lateral support extension	10271	71,00
<input type="radio"/> Small (14x9)		05642 188,00
<input type="radio"/> Large (18x12)		05646 188,00

- Headrest adjustable - Double ball mounting system
Available only with backrest Physio Posture





<input type="checkbox"/> Large (34x14)		05648 295,00
<input type="checkbox"/> Small (28x14)		

Date	Order	User	Address
------	-------	------	---------

Note

Wheelchair : Progeo TEKNA ADVANCE - Rev. 2024.1	Prices - VAT excluded <small>to be added to the price list</small>	<p style="text-align: center; margin: 0;">Send this page only if any optional have been chosen</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10px; text-align: center;">☆</td> <td style="width: 70%;">standard</td> <td style="width: 20%;"><small>(supplied if no requests are indicated)</small></td> </tr> <tr> <td style="text-align: center;">□</td> <td>on request</td> <td><small>(without extra charge)</small></td> </tr> <tr> <td style="text-align: center;">○</td> <td>optional</td> <td><small>(cost to be added at the basic price)</small></td> </tr> </table>	☆	standard	<small>(supplied if no requests are indicated)</small>	□	on request	<small>(without extra charge)</small>	○	optional	<small>(cost to be added at the basic price)</small>
☆	standard	<small>(supplied if no requests are indicated)</small>									
□	on request	<small>(without extra charge)</small>									
○	optional	<small>(cost to be added at the basic price)</small>									

REFERENCE TABLES

1	Front height			
	With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster			
	Casters			
FORK	3"	4"	5"	6"
SMALL 	max 49	max 50	max 51	
MEDIUM 	max 53	max 54	max 55	max 57
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
	Rear wheels			
	22"	24"	25"	26"
H min	min. 36	min. 38	min. 39	min. 40
H max	max 44	max 48	max 49	max 51

3	Footplate option	Footplate distance
	Automatic Closure footplate	min. 36
	Automatic Closure footplate , with multifunction rotative system	min. 31
	One piece footplate - Flip-up!	min. 36
	One piece footplate - Flip-up! with multifunction rotative system	min. 31
	One piece aluminium footplate - Flip-up!	min. 36
	Plastic footrest	min. 36
	One piece Carbon fibre footplate	min. 36
Calculated with the best configuration and with the footplates positioned on inside.		
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)