

Tilty Vario Tilttable wheelchair adolescents and adults



Tilty Vario



Tilttable wheelchair for
adolescents and adults
Active, tilttable and growable




Unique people, unique solutions.

Tilty Vario

Active, tiltable, growable


German aid index N° 18.99.02.1013



seat width
28-44 cm
(+2 cm)



back height
30-50 cm
(+5 cm)



seat depth
30-44 cm
(+4 cm)



weight capacity
100 kg



wheel sizes
22" and 24"



tilting angle
from **0° to +25°**



Overview

- Tilttable wheelchair for adolescents and adults
- Back adjustable in angle (80° – 120°) and completely stowable
- Highly stable rigid frame with 10 cm abduction; also suitable for Pörnbacher therapy methods
- Smooth handling with minimized effort
- Suitable as a seat for transport in a vehicle according to ISO 7176-19



www.sorgrollstuhltechnik.de

Equipment examples

All equipment options are on the order form of the Tilty Vario



holder to attach
abduction wedge




truss pads



arm pads adjustable in
height and angle



undivided foot rest



head rests with chest-
shoulder-belt

And
much
more ...