

## FOLDABLE

Public price starting from £ 3,700.00

M.D.R.: **2544458**

Basid UDI-DI: **805378246PGH02IBX**

Order  Date

Quote

Dealer

Order reference

## ALWAYS AT HAND



### MODEL FEATURES

- Triride Software
- Customisable electronics
- Display with 5 levels of power and speed
- Electronic brake
- Reverse
- Cruise Control
- Assisted load
- Weight: 8.2kg
- Range: up to 35 km
- Torque: 40Nm
- Lithium-ion cylindrical battery
- Tyre City or Cross
- Wheel: 12"
- ECODRIVE-Energy saving mode
- Charging the battery during downhill and braking

Small, light and compact, thanks to the sophisticated Made in Triride electronics, it can tackle various routes and is perfect indoors and for those who love to travel independently.

It allows a freedom of movement that is unmatched with other aids and once folded, it can be transported with the dedicated trolley.

### Notes


### WHEELCHAIR FEATURES

Brand

Model

Front tube's diameter

Seat width

Seat depth (E)

Toe to front tube distance (A)

Ground/knee height (D)

Length of stretched arms (C)

Front tubes angle ( degree )

Frame length (B)

Front tubes width(F)\*

User's weight

Footrest

Fixed  Internal

Removable  External

### USER'S PERSONAL DATA

Name

Tel.:

E-mail:

Age  Pathology:

### DELIVERY ADDRESS

Delivery to:

Postcode


Address

Town

County

\* The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm

**All fields are obligatory**



**N.B.:** For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc )

Assembly, installation and instruction for use carried out by a Triride specialist technician  (on request)

## MTW: TREKKING POWER

	Adaptation for MTW Trekking Power Kit purchase ( complete and attach dedicated form )	no cost	<input type="checkbox"/>
T0013	Adaptation for MTW Trekking Power ( future purchase )	£185.00	<input type="checkbox"/>

## TRIRIDE CONFIGURATION

### Electronics

BS-O-096	Performance Pack : IBS+ICC+IBS button	£171.00	<input type="checkbox"/>
BS-O-101	IBS - Intelligent Braking System	£68.00	<input type="checkbox"/>
BS-O-097	ICC - Intelligent Cruise Control	£68.00	<input type="checkbox"/>
BS-O-086	ECODRIVE - Increased battery range	£68.00	<input type="checkbox"/>
BS-O-113	MDC - Motion Direct Control	£132.00	<input type="checkbox"/>
BS-R-366	Electronic brake button / IBS	£44.00	<input type="checkbox"/>
BS-O-076	Analogue brake with thumb control	£95.00	<input type="checkbox"/>
BS-O-077	Automatic brake on accelerator release	£41.00	<input type="checkbox"/>
BS-O-079	Emergency brake with remote control	£95.00	<input type="checkbox"/>

### Display LCD

BS-C-076	Display LCD Basic	default	<input type="checkbox"/>
BS-C-077	Display LCD Special with USB	£35.00	<input type="checkbox"/>
BS-C-078	Display LCD Colour with USB	£86.00	<input type="checkbox"/>
BS-C-079	Display Colour with USB with customisable Triride programming	£142.00	<input type="checkbox"/>
	customisation: <input type="text"/>	( 21 character )	

### Tyre

BS-R-353	Tyre 12" x2.125 city long life	£0.00	<input type="checkbox"/>
BS-R-297	Tyre 12" x2.125 cross long life	£11.00	<input type="checkbox"/>

### Throttle

BS-R-299	Thumb throttle ambidextrous	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£0.00
L: BS-R-349 R: BS-R-348	Palm throttle with long lever	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£20.00
L: BS-R-333 R: BS-R-335	Half twist throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£38.00
L: T0062 R: BS-R-324	Full twist throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£38.00

### Brakes

L: BS-R-319 R: BS-R-318	Single mechanical brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£0.00
L: BS-O-203 R: BS-O-273	Single mechanical parking brake (not available with hydraulic brake)	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£14.00
L: BS-O-274 R: BS-O-275	Single hydraulic brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£122.00
BS-O-218	Mechanical brake - double caliper		£149.00
BS-O-194	Hydraulic brake - double caliper		£224.00

N.B: The electronic brake ( including IBS when required) can be actioned also by the brake levers of the device

### Frame colour

BS-O-078	Titanium	default
BS-O-078-P	Polished	£252.00
BS-O-217-N	Black Ral 9005	£135.00
BS-O-217-R	Red Ral 3020	£135.00
BS-O-217-W	White Ral 9003	£135.00
BS-O-217-G	Yellow Ral 1003	£135.00
BS-O-217-A	Orange Ral 2004	£135.00
BS-O-217-B	Blue Ral 5017	£135.00
BS-O-217	RAL on request	£207.00

### Personalized frame wraps

BS-O-087-M	Camouflage	£504.00
BS-O-087-C	Carbon Effect	£504.00
BS-O-087-F	Graffiti wrap Black/White	£504.00
	Personalisation: Two-tone, gradients, etc.	on request



**Mudguard**

BS-O-245	Mudguard Black	default	<input type="checkbox"/>
BS-O-167	Mudguard Blue	£0.00	<input type="checkbox"/>
BS-O-169	Mudguard Orange	£0.00	<input type="checkbox"/>
BS-O-166	Mudguard Green	£0.00	<input type="checkbox"/>
BS-O-164	Mudguard Red	£0.00	<input type="checkbox"/>
BS-O-215	Mudguard Carbon	£96.00	<input type="checkbox"/>

**Handlescrews**

BS-O-301	Handlescrews - Black	£0.00	<input type="checkbox"/>
BS-O-128	Handlescrews - Blue	£0.00	<input type="checkbox"/>
BS-O-127	Handlescrews - Orange	£0.00	<input type="checkbox"/>
BS-O-130	Handlescrews - Green	£0.00	<input type="checkbox"/>
BS-O-129	Handlescrews - Red	£0.00	<input type="checkbox"/>
BS-O-139	Lever handlescrews for Tetraplegic user	£5.00	<input type="checkbox"/>

**SPEED SETTING**

We herewith advise, that using the Triride with max. speed settings higher than 4mph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >4mph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 4mph.

**When placing an order for a max. speed-setting higher than 4mph please confirm with your signature below**, that you have understood the above information, you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Up to 4mph (6 Kph)

Up to 8mph (13 Kph)

>8mph (13kph)

Date

Customer Name

Customer signature

**OPTIONAL**

**Attachment Kit for SECOND WHEELCHAIR - please provide details below ( prices vary as per specifics )**

Tilting frame

Connecting arms

Clamps

Brand

Model

Seat width  cm

Seat depth  cm

Front tubes diameter  mm

## ACCESSORIES / OPTIONAL

BS-R-344	Battery 36V 7A (252Wh) Cylindrical		£468.00	
BS-O-012	Neoprene 20/19 battery cover		£24.00	
BS-O-083	Central arch stand		£149.00	
BS-O-278	Lateral stand with wheels		£113.00	
BS-O-229	Rearview mirror	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£29.00	
BS-O-251	Carry Trailer for fixed frame	Rear axle diameter <input type="text"/> mm	£279.00	
BS-O-214	Leg cover		£19.00	
BS-O-292	Calf support with adaptation to tilting frame		£142.00	
BS-O-115	Triride bag with Trolley		£162.00	
BS-O-114	Arms stopper rings, pair		£16.00	
BS-O-174	Laser customised engraving on frame (on inox and black frame only)		£68.00	
BS-O-156	Digital engraving		£68.00	
BS-O-157	Personalised sticker		£59.00	
BS-O-081	Integrated front light, activated from display		£95.00	
BS-O-021	Front and rear light kit with switch on handlebar		£162.00	
BS-O-232	Light - Front 1200 Lumen		£100.00	
BS-O-198	Light Front USB		£29.00	
BS-O-180	Light - Rear USB		£33.00	
BS-O-015	Rear axle extension for wheelchair ( please complete dedicated form )		on request	

Prices are intended for public and before tax

**For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website [web www.triride.uk.com](http://www.triride.uk.com)**

**N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product**