

# **KIDS**

#### Public price starting from £ 3,495.00

M.D.R.: **2544454** 

Basid UDI-DI: **805378246KDS014BD** 

Order

Date

Quote

Dealer

Order reference

## **MODEL FEATURES**

- · Triride Software
- · Customisable electronics
- · Display with 5 levels of power and speed
- · Electronic brake
- · Reverse
- · Cruise Control
- · Assisted load
- · Weight: 7.5kg
- · Range: up to 35 km
- · Torque: 40Nm
- · Lithium-ion cylindrical battery
- · Tyre City or Cross
- · Wheel: 10"or 12"
- · ECODRIVE-Energy saving mode
- $\boldsymbol{\cdot}$  Charging the battery during downhill and braking

## FOR THE LITTLE CHAMPIONS!




Triride for the little ones, derived from the Light model. Size adjustable for developmental age. The 10" or 12" wheel makes it super compact, very light and easy to handle, ideal for both indoor and outdoor spaces. Power and performance controlled by Triride Electronics.

Notes		



WHEELCHAIR FEATURES		USER'S PERSONAL DATA			
Brand		Name			
Model		Tel:			
Front tube's diameter		E-mail:			
Seat width		Age Pathology:			
Seat depth (E)					
Toe to front tube distance (A)		DELIVERY ADDRESS			
Ground/knee height (D)		DELIVERY ADDRESS			
Length of stretched arms (C)		Delivery to:			
Front tubes angle ( degree )					
Frame length (B)		Postcode			
Front tubes width(F)*		Address			
User's weight					
Footrest		Town			
Fixed Internal		County			
Removable	External				

<sup>\*</sup> The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm All fields are obligatory



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc.)

(on request) Assembly, installation and instruction for use carried out by a Triride specialist technician



	MTW: TREKKING POWER	
Adaptation for M	TW Trekking Power Kit purchase ( complete and attach dedicated form )	no cost
T0013	Adaptation for MTW Trekking Power (future purchase)	£185.00
	TRIRIDE CONFIGURATION	
Electronics		
BS-O-096	Performance Pack : IBS+ICC+IBS button	included
BS-O-113	MDC - Motion Direct Control	£132.00
BS-O-076	Analogue brake with thumb control	£95.00
BS-O-077	Automatic brake on accelerator release	£41.00
BS-O-079	Emergency brake with remote control	£95.00
Display LCD		
BS-C-076	Display LCD Basic	default
BS-C-077	Display LCD Special with USB	£35.00
BS-C-078	Display LCD Colour with USB	£86.00
BS-C-079	Display Colour with USB with customisable Triride programming	£142.00
	customisation: (2)	l character )
Tyre		_
BS-R-280	Wheel 10" with 10"x2.5 city tyre	£162.00
BS-R-353	Tyre 12" x2.125 city long life	£0.00
BS-R-297	Tyre 12" x2.125 cross long life	£0.00



Throttle				
BS-R-299	Thumb throttle ambidextrous	LHS	RHS	£0.00
L: BS-R-349				
R: BS-R-348	Palm throttle with long lever	LHS	RHS	£20.00
L: BS-R-333		_	_	
R: BS-R-335	Half twist throttle	LHS	RHS	£38.00
Brakes				
L: BS-R-319			DUG	60.00
R: BS-R-318	Single mechanical brake	LHS	RHS	£0.00
L: BS-O-203 R: BS-O-273	Single mechanical parking brake	LHS	RHS	£14.00
	(not available with hydraulic brake)	LIIS	CIIN	L14.00
L: BS-O-274 R: BS-O-275	Single hydraulic brake (not for 10")	LHS	RHS	£122.00
	: brake ( including IBS when required) can be actioned also by t			
Handlebar	, arake (melaanig 125 when required) can be accioned also by c	ine brane n	evers of the device	
BS-O-246	Straight (indicate width between 46cm and 50cr	m)	cm	£0.00
BS-O-208	4 Bends riser	£0.00		
BS-O-107	4 Bends high riser	£0.00		
Handlebar sto				
BS-O-222	Fixed handlebar stem			£0.00
BS-O-293	Collapsible handlebar	£44.00		
Custom frame	9			
BS-C-817	Height adjustable telescopic stem			£108.00
Frame colour				
BS-O-078	Titanium			default
BS-O-078-P	Polished			£252.00
BS-O-217-N	Black Ral 9005			£0.00
BS-O-217-R	Red Ral 3020			£0.00
BS-O-217-W	White Ral 9003			£0.00
	Yellow Ral 1003			£0.00
BS-O-217-G				
	Orange Ral 2004			£0.00
BS-O-217-G BS-O-217-A BS-O-217-B	Orange Ral 2004 Blue Ral 5017			£0.00



Persona	lized fra	ame wraps
---------	-----------	-----------

BS-O-087-M	Camouflage	£504.00
BS-O-087-C	Carbon Effect	£504.00
BS-O-087-F	Graffiti wrap Black/White	£504.00
	Personalisation: Two-tone, gradients, etc.	on request
Mudguard		
BS-O-245	Mudguard Black	default
BS-O-167	Mudguard Blue	£0.00
BS-O-169	Mudguard Orange	£0.00
BS-O-166	Mudguard Green	£0.00
BS-O-164	Mudguard Red	£0.00
BS-O-215	Mudguard Carbon	£96.00
Handlescrews	S	
BS-O-301	Handlescrews - Black	£0.00
BS-O-128	Handlescrews - Blue	£0.00
BS-O-127	Handlescrews - Orange	£0.00
BS-O-130	Handlescrews - Green	£0.00
BS-O-129	Handlescrews - Red	£0.00
BS-O-139	Lever handlescrews for Tetraplegic user	£5.00

#### **SPEED SETTING**

We herewith advise, that using the Triride with max. speed settings higher than 4mph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >4mph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 4mph.

When placing an order for a max. speed-setting higher than 4mph please confirm with your signature below, that you have understood the above information, you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Up to 4mp	h (6 Kph )	Up to 8m	nph (13 Kph)	>	8mph (13kph)	
Date		Customer Name			Customer signa	ature
			OPTIONA	<b>NL</b>		
Attachment Kit fo	Attachment Kit for SECOND WHEELCHAIR - please provide details below (prices vary as per specifics)					cifics )
Tilting frame		Connecting arms		Clamps		
Brand		Model				
Seat width	cm	Seat depth	cm	Front tubes diamete	er mr	n



ACCESSORIES / OPTIONAL				
BS-R-344	Battery 36V 7A (252Wh) Cylindrical	£468.00		
BS-O-012	Neoprene 20/19 battery cover	£24.00		
BS-O-083	Central arch stand	£149.00		
BS-O-278	Lateral stand with wheels	£113.00		
BS-O-229	Rearview mirror LHS RHS	£29.00		
BS-O-251	Carry Trailer for fixed frame Rear axle diametermm	£279.00		
BS-O-214	Leg cover	£19.00		
BS-O-292	Calf support with adaptation to tilting frame	£142.00		
BS-O-115	Triride bag with Trolley	£162.00		
BS-O-114	Arms stopper rings, pair	£16.00		
BS-O-174	Laser customised engraving on frame (on inox and black frame only)	£68.00		
BS-O-156	Digital engraving	£68.00		
BS-O-157	Personalised sticker	£59.00		
BS-O-081	Integrated front light, activated from display	£95.00		
BS-O-021	Front and rear light kit with switch on handlebar	£162.00		
BS-O-232	Light - Front 1200 Lumen	£100.00		
BS-O-198	Light Front USB	£29.00		
BS-O-180	Light - Rear USB	£33.00		
BS-O-015	Rear axle extension for wheelchair ( please complete dedicated form )	on request		

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website web www.triride.uk.com

N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product