

# Triride®

by Gianni Conte

A NEW FREEDOM



USB

## mod. SPECIAL L14



CE



Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC



14"



36 V / 48 V



1200 W / 1500 W



12,4 Kg<sup>(1)</sup>



up to 50 Km<sup>(2)</sup>

### TECHNICAL FEATURES

EN

Max Power 36V 1200W / 48V 1500W

Electronics by Triride

Reverse

Cruise Control

Assisted Load

Electronic Brake

Weight: 12,4 Kg

14" driving wheel

LCD Display

36V or 48V Lithium-ion battery

Powerful disc brake

Tire selection: City or Cross

Stainless steel frame

Speed compliant with legal standards

5 Speed levels

Customized frame color

<sup>1</sup> Weight refers to device excluding battery and arms for wheelchair connection.

<sup>2</sup> Autonomy refers to device with EcoDrive installed, 48 V battery on a flat route at medium speed of 6 km/h (level 2) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.