

MODEL SPECIAL LIGHT

Medical Devices Registration N: **1560509**

ORDER

QUOTE

ORDER REFERENCE

DATE

Please complete and send this form by e-mail to: orders@triride.uk

SPECIAL LIGHT **£4,150.00**



12" City tyre long life £27.00
 12" Cross tyre long life £36.00
 14" City tyre long Life standard

SPECIAL LIGHT ALLROAD **£4,300.00**



With 3.00 - 8 ALL ROAD tyre

CODE
TR-SP-LIGHT


DESCRIPTION
Wheel 12" or 14", Programmable and customizable Triride electronics, reverse, cruise control, assisted lifting, dual map: normal/sport, 5 speeds, electronic brake, display LCD, battery recharging during braking

USER'S PERSONAL DATA****

Name _____
 Surname _____
 Phone _____
 E-mail _____
 Age _____
 Pathology _____

DELIVERY ADDRESS

Delivery to _____
 Street _____
 Town _____
 Postcode _____

WHEELCHAIR FEATURES																					
Brand	_____																				
Model	_____																				
Front tubes diameter	_____ mm																				
Seat width	_____ cm																				
Seat depth (b)	_____ cm																				
Toe to front tube distance (a)*	_____ cm																				
Ground/knee height (d)	_____ cm																				
Length of arms from the shoulder back (c)	_____ cm																				
Front tubes angle	_____ degrees																				
User's weight	_____ kg																				
Footrest	Fixed <input type="checkbox"/> Internal <input type="checkbox"/> Removable <input type="checkbox"/> External <input type="checkbox"/>																				
For wheelchairs with removable footrests, measure as per "a" from toe to front tube.																					
																					
TRIRIDE CONFIGURATION																					
Electronics	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">BS-O-096</td> <td style="width: 55%;">Performance KIT=Autonomy and Security(IFS+ACC+Button)</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 20%; text-align: right;">£205.00</td> </tr> <tr> <td>BS-O-101</td> <td>IBS-Intelligent Braking System (included with push&pull)</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£82.00</td> </tr> <tr> <td>BS-O-097</td> <td>ACC-Adaptive Cruise Control (included with push&pull)</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£82.00</td> </tr> <tr> <td>BS-R-366</td> <td>Electronic brake button/IBS <input type="checkbox"/> LHS <input type="checkbox"/> RHS</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£100.00</td> </tr> <tr> <td>BS-O-113</td> <td>MDC- Motion Direct Control</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£132.00</td> </tr> </table>	BS-O-096	Performance KIT=Autonomy and Security(IFS+ACC+Button)	<input type="checkbox"/>	£205.00	BS-O-101	IBS-Intelligent Braking System (included with push&pull)	<input type="checkbox"/>	£82.00	BS-O-097	ACC-Adaptive Cruise Control (included with push&pull)	<input type="checkbox"/>	£82.00	BS-R-366	Electronic brake button/IBS <input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£100.00	BS-O-113	MDC- Motion Direct Control	<input type="checkbox"/>	£132.00
BS-O-096	Performance KIT=Autonomy and Security(IFS+ACC+Button)	<input type="checkbox"/>	£205.00																		
BS-O-101	IBS-Intelligent Braking System (included with push&pull)	<input type="checkbox"/>	£82.00																		
BS-O-097	ACC-Adaptive Cruise Control (included with push&pull)	<input type="checkbox"/>	£82.00																		
BS-R-366	Electronic brake button/IBS <input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£100.00																		
BS-O-113	MDC- Motion Direct Control	<input type="checkbox"/>	£132.00																		
Special commands for tetraplegic users	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">BS-O-175</td> <td style="width: 55%;">Push & Pull for Tetraplegic users (brake/accelerator push & pulls)</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/> LHS <input type="checkbox"/> RHS</td> <td style="width: 20%; text-align: center;"><input type="checkbox"/></td> <td style="width: 10%; text-align: right;">£387.00</td> </tr> <tr> <td>BS-O-176</td> <td>Push & Pull for Tetraplegic users double brake/accelerator push & pull</td> <td style="text-align: center;"><input type="checkbox"/> LHS <input type="checkbox"/> RHS</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£441.00</td> </tr> <tr> <td>BS-O-307</td> <td>Fixed vertical extension</td> <td></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£15.00</td> </tr> </table> <p style="text-align: right;">Configuration N. _____</p> <p>"For more information please see the here attached page with different configurations"</p>	BS-O-175	Push & Pull for Tetraplegic users (brake/accelerator push & pulls)	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£387.00	BS-O-176	Push & Pull for Tetraplegic users double brake/accelerator push & pull	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£441.00	BS-O-307	Fixed vertical extension		<input type="checkbox"/>	£15.00					
BS-O-175	Push & Pull for Tetraplegic users (brake/accelerator push & pulls)	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£387.00																	
BS-O-176	Push & Pull for Tetraplegic users double brake/accelerator push & pull	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/>	£441.00																	
BS-O-307	Fixed vertical extension		<input type="checkbox"/>	£15.00																	
Inner tube	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">BS-R-300</td> <td style="width: 55%;">12" anti puncture (flatfree) not available for ALL ROAD tyre</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 20%; text-align: right;">£9.00</td> </tr> <tr> <td>BS-R-342</td> <td>14" anti puncture (flatfree) not available for ALL ROAD tyre</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£11.00</td> </tr> </table>	BS-R-300	12" anti puncture (flatfree) not available for ALL ROAD tyre	<input type="checkbox"/>	£9.00	BS-R-342	14" anti puncture (flatfree) not available for ALL ROAD tyre	<input type="checkbox"/>	£11.00												
BS-R-300	12" anti puncture (flatfree) not available for ALL ROAD tyre	<input type="checkbox"/>	£9.00																		
BS-R-342	14" anti puncture (flatfree) not available for ALL ROAD tyre	<input type="checkbox"/>	£11.00																		
Display LCD	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">BS-C-077</td> <td style="width: 55%;">Display LCD Special with USB</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 20%; text-align: right;">standard</td> </tr> <tr> <td>BS-C-078</td> <td>Display LCD Color USB</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: right;">£50.00</td> </tr> </table>	BS-C-077	Display LCD Special with USB	<input type="checkbox"/>	standard	BS-C-078	Display LCD Color USB	<input type="checkbox"/>	£50.00												
BS-C-077	Display LCD Special with USB	<input type="checkbox"/>	standard																		
BS-C-078	Display LCD Color USB	<input type="checkbox"/>	£50.00																		

Battery	BS-R-369	Battery 36V 7A (252 Wh)	<input type="checkbox"/>	standard
	BS-R-394	Battery 36V 8.8 A (316Wh)	<input type="checkbox"/>	standard
	BS-R-393-S	Battery 36V 14 A (504 Wh)	<input type="checkbox"/>	£155.00
	BS-R-368-S	Battery 48V 5.2 A (249,6 Wh)	<input type="checkbox"/>	£110.00
	BS-R-395-S	Battery 48V 10.5 A (504 Wh)	<input type="checkbox"/>	£200.00
Handlebar	BS-O-246	Straight	<input type="checkbox"/>	standard
	BS-O-208	4 bends Zoom	<input type="checkbox"/>	standard
	BS-O-107	4 bends	<input type="checkbox"/>	standard
	BS-O-207	Special handlebar for tetra users	<input type="checkbox"/>	£152.00
Handlebar stem	BS-O-222	Adjustable 110 mm	<input type="checkbox"/>	standard
	BS-O-220	Adjustable 130 mm	<input type="checkbox"/>	standard
	BS-O-185	Adjustable 160 mm	<input type="checkbox"/>	£76.00
	BS-O-211	Adjustable 190 mm	<input type="checkbox"/>	£76.00
Throttle	BS-R-299	Thumb throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/> standard
	LS:S-R-329 R-328	RS:S- Palm throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/> standard
	S-R-333	Half twist throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/> £36.00
Brakes	LS:S-R-331 R-306	RS:S- Single mechanical disk brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/> standard
	BS-O-218	Mechanical disk brake - double caliper		<input type="checkbox"/> £135.00
	LS:BS-O-274 RS:BS-O-275	Single hydraulic disk brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	<input type="checkbox"/> £117.00
	BS-O-194	Hydraulic disk brake - double caliper		<input type="checkbox"/> £269.00
N.B.:The electronic brake (including IBS when required) can be actioned also only by the brake levers of the Triride.				
Frame colour	BS-C-078	Neutral (titanium steel colour)	<input type="checkbox"/>	standard
	BS-C-078-P	Polished	<input type="checkbox"/>	£177.00
	BS-O-217-N	Black Ral 9005	<input type="checkbox"/>	£132.00
	BS-O-217-R	Red Ral 3020	<input type="checkbox"/>	£132.00
	BS-O-217-W	White Ral 9003	<input type="checkbox"/>	£132.00
	BS-O-217-G	Yellow Ral 1023	<input type="checkbox"/>	£132.00
	BS-O-217-A	Orange Ral 2004	<input type="checkbox"/>	£132.00
	BS-O-217-B	Blu Ral 5015	<input type="checkbox"/>	£132.00
	BS-O-217	RAL on request RAL N. _____	<input type="checkbox"/>	£177.00
	BS-O-087	Custom graphic design upon request	<input type="checkbox"/>	
	BS-O-087-M	Camouflage	<input type="checkbox"/>	£502.00
	BS-O-087-C	Carbon effect	<input type="checkbox"/>	£502.00
	BS-O-087-F	Graffiti graphics	<input type="checkbox"/>	£502.00

Mudguard	BS-O-245	PP Black(not available on Allroad)	<input type="checkbox"/>	standard
	BS-O-167	PP Blue (not available for Allroad)	<input type="checkbox"/>	standard
	BS-O-169	PP Orange (not available for Allroad)	<input type="checkbox"/>	standard
	BS-O-166	PP Green (not available for Allroad)	<input type="checkbox"/>	standard
	BS-O-164	PP Red (not available for Allroad)	<input type="checkbox"/>	standard
	BS-O-215	Composite carbon (standard on Allroad)	<input type="checkbox"/>	£94.00
Triangular handlescrews	BS-O-301	Triangular handlescrews - black	<input type="checkbox"/>	standard
	BS-O-130	Triangular handlescrews - green	<input type="checkbox"/>	standard
	BS-O-128	Triangular handlescrews - blue	<input type="checkbox"/>	standard
	BS-O-127	Triangular handlescrews - orange	<input type="checkbox"/>	standard
	BS-O-129	Triangular handlescrews - red	<input type="checkbox"/>	standard
Lever handlescrews	BS-R-330	Lever handlescrews	<input type="checkbox"/>	standard

OPTIONAL

SECOND WHEELCHAIR FEATURES FOR EXTRA CLAMPS			NOTES			
Tilting frame+arms	<input type="checkbox"/>		<hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/> <hr/>			
Brand _____	Model _____	_____				
Front tubes diameter		_____ mm				
Seat width		_____ cm				
Seat depth (b)		_____ cm				
Toe to front tube distance (a)		_____ cm				
Ground/knee hight (d)		_____ cm				
Lenght of arms from the shoulder back (c)		_____ cm				
Front tubes angle		_____ degrees				
Footrest	Fixed	<input type="checkbox"/>			Internal	<input type="checkbox"/>
	Removable	<input type="checkbox"/>			Esternal	<input type="checkbox"/>
* For wheelchairs with removable footrests, measure as per "a" from toe to front tube.						

ACCESSORIES				
36V 7A Battery (252 Wh)			BS-R-369	<input type="checkbox"/> £544.00
36V 8.8 A Battery (316Wh)			BS-R-394	<input type="checkbox"/> £544.00
36V 14 A Battery (504 Wh)			BS-R-393	<input type="checkbox"/> £700.00
48V 5.2 A Battery (249,6 Wh)			BS-R-368	<input type="checkbox"/> £650.00
48V 10.5 A Battery (504 Wh)			BS-R-395	<input type="checkbox"/> £753.00
Neoprene battery cover			BS-O-012	<input type="checkbox"/> £23.00
Front spare battery holder			BS-R-294	<input type="checkbox"/> £157.00
T-Stand with wheels			BS-O-201	<input type="checkbox"/> £109.00
Rear mirror	<input type="checkbox"/> LHS	<input type="checkbox"/> RHS	BS-O-229	<input type="checkbox"/> £26.00
**Carry Trailer		Rear axle diameter _____mm	BS-O-251	<input type="checkbox"/> £269.00
Handlebar basket with fabric cover			BS-O-011	<input type="checkbox"/> £93.00
Leg cover			BS-O-214	<input type="checkbox"/> £15.00
Triride bag			BS-O-213	<input type="checkbox"/> £141.00
Arms stopper rings, pair			BS-O-114	<input type="checkbox"/> £14.00
Laser customized engraving on the frame			BS-O-174	<input type="checkbox"/> £68.00
Light - Front 1200 Lumen			BS-O-232	<input type="checkbox"/> £100.00
Light - Front+USB			BS-O-198	<input type="checkbox"/> £27.00
Light - Rear+USB			BS-O-180	<input type="checkbox"/> £32.00
"For a complete list of accessories and spare parts visit our website: www.triride.uk "				
* Mandatory fields for warranty purposes				
**Available for fixed frame wheelchairs, communicate the diameter of the rear axle of the wheelchair				
***Prices are intended for public and before tax				
N. B.:The Company reserves the exclusive right to change, without notice, if it deems is useful for the improvement of the Product				

DISPLAY LCD		
		
BS-C-076	BS-C-077	BS-C-078
HANDLEBAR		
		
BS-O-246	BS-O-208	BS-O-107
		
BS-O-099		
ACCELERATOR		
		
S-R-328/S-R-329	BS-R-299	S-R-333
DOUBLE BATTERY SUPPORT	LONG LIFE TIRE	
		
BS-R-294	CITY	CROSS

OPTIONAL



BS-O-201



BS-O-229



BS-O-251



BS-O-114



BS-O-214



BS-O-213



FENDER PP



BS-O-215



BS-O-232



BS-O-198



BS-O-180



BS-O-012

OPTIONAL



KNOB SCREWS



BS-R-330



~~BS-O-11~~

Obsolete



BS-O-011

Basket with fabric cover
3 external pockets
1 inside pocket
Rain cover



BS-O-174

CUSTOMIZATION



POLISH BS-C-078-P



CUSTOM GRAPHICS BS-O-087



Simil Carbon



Comics