

T-ROCKS

Public price starting from £5,995.00

MDR: **2544507**

Basic UDI-DI: **805378246TKSO76J5**

Order Date

Quote

Dealer

Order reference

MODEL FEATURES

- Triride Software
- Customisable electronics
- Display with 5 levels of power and speed
- Electronic brake
- Reverse
- Cruise Control
- Assisted load
- Weight: 12.8 kg
- Range: up to 50 km
- Torque: 65Nm
- Lithium-ion battery 48V 10.5A
- Wheel size 16" with Carbon mudguard
- Tyre 16"x4 cross
- ECODRIVE-Energy saving mode
- Two driving modes: Normal/Sport
- Charging the battery during downhill and braking

Go beyond all limits with T-Rocks. Designed to accompany you on the most demanding outdoor excursions in contact with nature. It is ideal for Off-Road routes and will give you unexpected emotions. T-Rocks evolves with a new fork design, even more performing and expands the family with a model with a 16"x4 wheel which combined with the exclusive HT motor guarantees excellent traction and power to tackle all terrains. Available with double battery (optional) with selection switch. IBS and ICC as standard which allow professional assistance and complete customization of performance.

FOR OFF-ROAD LOVERS



WHEELCHAIR FEATURES

Brand

Model

Front tube's diameter

Seat width

Seat depth (E)

Toe to front tube distance (A)

Ground/knee height (D)

Length of stretched arms (C)

Front tubes angle (degree)

Frame length (B)

Front tubes width(F)*

User's weight

Footrest

Fixed Internal

Removable External

USER'S PERSONAL DATA

Name

Tel.:

E-mail:

Age Pathology:

DELIVERY ADDRESS

Delivery to:

Postcode


Address

Town

County

* The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm

All fields are obligatory



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc)

Assembly, installation and instruction for use carried out by a Triride specialist technician (on request)



Notes

MTW: TREKKING POWER

Adaptation for MTW Trekking Power Kit purchase (complete and attach dedicated form)	no cost	<input type="checkbox"/>
T0013 Adaptation for MTW Trekking Power (future purchase)	£185.00	<input type="checkbox"/>

COMMANDS FOR TETRAPLEGIC USERS

Configuration for Tetraplegic Commands N. <input type="text"/>	<input type="checkbox"/>
(For the desired configuration fill and attach the dedicated form)	<input type="checkbox"/>
Personalised configuration (attach picture or drawing)	<input type="checkbox"/>
"For more information, see the attached page with various configurations"	

TRIRIDE CONFIGURATION

Battery

BS-R-368-S Battery 48V 6.2 USB (297 Wh)	£0.00	<input type="checkbox"/>
BS-R-395-S Battery 48V 10.5 USB (504 Wh)	default	<input type="checkbox"/>

Electronics

BS-O-096 Performance Pack : IBS+ICC+IBS button	included	<input type="checkbox"/>
BS-O-113 MDC - Motion Direct Control	£132.00	<input type="checkbox"/>
BS-O-076 Analogue brake with thumb control	£95.00	<input type="checkbox"/>
BS-O-077 Automatic brake on accelerator release	£41.00	<input type="checkbox"/>

Display LCD

BS-C-077 Display LCD Special with USB	default	<input type="checkbox"/>
BS-C-078 Display LCD Colour with USB	£53.00	<input type="checkbox"/>
BS-C-079 Display Colour with USB with customisable Triride programming	£126.00	<input type="checkbox"/>
customisation: <input type="text"/>	(21 character)	



Throttle

BS-R-299	Thumb throttle ambidextrous	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£0.00	<input type="checkbox"/>
L: BS-R-349				
R: BS-R-348	Palm throttle with long lever	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£20.00	<input type="checkbox"/>
L: BS-R-333				
R: BS-R-335	Half twist throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£38.00	<input type="checkbox"/>
L: T0062				
R: BS-R-324	Full twist throttle	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£38.00	<input type="checkbox"/>

Brakes

L: BS-R-319				
R: BS-R-318	Single mechanical brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£0.00	<input type="checkbox"/>
L: BS-O-203				
R: BS-O-273	Single mechanical parking brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£14.00	<input type="checkbox"/>
L: BS-O-274	(not available with hydraulic brake)			
R: BS-O-275	Single hydraulic brake	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£122.00	<input type="checkbox"/>
BS-O-218	Mechanical brake - double caliper		£149.00	<input type="checkbox"/>
BS-O-194	Hydraulic brake - double caliper		£224.00	<input type="checkbox"/>

N.B: The electronic brake (including IBS when required) can be actioned also by the brake levers of the device

Handlebar

BS-O-099	Over MTB		default	<input type="checkbox"/>
----------	----------	--	---------	--------------------------

Handlebar stem

BS-O-098	Adjustable 120mm stem		default	<input type="checkbox"/>
----------	-----------------------	--	---------	--------------------------

Frame colour

BS-O-078	Titanium		default	<input type="checkbox"/>
BS-O-078-P	Polished		£252.00	<input type="checkbox"/>
BS-O-217-N	Black Ral 9005		£135.00	<input type="checkbox"/>
BS-O-217-R	Red Ral 3020		£135.00	<input type="checkbox"/>
BS-O-217-W	White Ral 9003		£135.00	<input type="checkbox"/>
BS-O-217-G	Yellow Ral 1003		£135.00	<input type="checkbox"/>
BS-O-217-A	Orange Ral 2004		£135.00	<input type="checkbox"/>
BS-O-217-B	Blue Ral 5017		£135.00	<input type="checkbox"/>
BS-O-217	RAL on request		£207.00	<input type="checkbox"/>



Personalized frame wraps

BS-O-087-M	Camouflage	£504.00	<input type="checkbox"/>
BS-O-087-C	Carbon Effect	£504.00	<input type="checkbox"/>
BS-O-087-F	Graffiti wrap Black/White	£504.00	<input type="checkbox"/>
	Personalisation: Two-tone, gradients, etc.	on request	<input type="checkbox"/>

Handlescrews

BS-O-301	Handlescrews - Black	£0.00	<input type="checkbox"/>
BS-O-128	Handlescrews - Blue	£0.00	<input type="checkbox"/>
BS-O-127	Handlescrews - Orange	£0.00	<input type="checkbox"/>
BS-O-130	Handlescrews - Green	£0.00	<input type="checkbox"/>
BS-O-129	Handlescrews - Red	£0.00	<input type="checkbox"/>
BS-O-139	Lever handlescrews for Tetraplegic user	£5.00	<input type="checkbox"/>

SPEED SETTING

We herewith advise, that using the Triride with max. speed settings higher than 4mph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >4mph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 4mph.

When placing an order for a max. speed-setting higher than 4mph please confirm with your signature below, that you have understood the above information, you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Up to 4mph (6 Kph)

Up to 8mph (13 Kph)

>8mph (13kph)

Date

Customer Name

Customer signature

OPTIONAL

Attachment Kit for SECOND WHEELCHAIR - please provide details below (prices vary as per specifics)

Tilting frame

Connecting arms

Clamps

Brand

Model

Seat width cm

Seat depth cm

Front tubes diameter mm

ACCESSORIES / OPTIONAL

BS-R-368	Battery 48V 6.2 USB (297 Wh)		£651.00	
BS-R-395	Battery 48V 10.5 USB (504 Wh)		£753.00	
BS-O-012	Neoprene 20/19 battery cover		£24.00	
BS-R-294	Front spare battery holder with switch		£224.00	
BS-O-083	Central arch stand		£149.00	
BS-O-278	Lateral stand with wheels		£113.00	
BS-O-229	Rearview mirror	<input type="checkbox"/> LHS <input type="checkbox"/> RHS	£29.00	
BS-O-251	Carry Trailer for fixed frame	Rear axle diameter <input type="text"/> mm	£279.00	
BS-O-019	Handlebar bag with B-Twin adaptor 10 L		£144.00	
BS-O-020	Handlebar bag with B-Twin adaptor 6 L		£144.00	
BS-O-011	Handlebar basket with fabric cover 24 L		£144.00	
BS-O-214	Leg cover		£19.00	
BS-O-292	Calf support with adaptation to tilting frame		£142.00	
BS-O-213	Triride bag		£153.00	
BS-O-114	Arms stopper rings, pair		£16.00	
BS-O-174	Laser customised engraving on frame (on inox and black frame only)		£68.00	
BS-O-156	Digital engraving		£68.00	
BS-O-157	Personalised sticker		£59.00	
BS-O-081	Integrated front light, activated from display		default	
BS-O-022	Double front spotlight, activated from display		£95.00	
BS-O-232	Light - Front 1200 Lumen		£100.00	
BS-O-198	Light Front USB		£29.00	
BS-O-180	Light - Rear USB		£33.00	
BS-O-015	Rear axle extension for wheelchair (please complete dedicated form)		on request	

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website [web www.triride.uk.com](http://www.triride.uk.com)

N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product