

TRIBIKE HYBRID PLUS

Public price starting from £ 4,550.00

Date

M.D.R: **2544521**

Basic UDI-DI: 805378246ETBKP113ZA

Order Quote

Dealer

Order reference

MODEL FEATURES

- · Triride Software
- · Performance Pack: IBS, ICC, throttle reverse, motor brake, regen battery
- ECODRIVE Energy saving mode
- · Display LCD with 5 speed levels
- · Electronic assisted start/climb
- · Battery USB 36V or 48V
- · Speed with motor support max 25km/h
- · Sensor assisted pedaling (PAS)
- · Weight: 11 kg
- ·Tyre: City or Cross
- · Wheel 16" or 20"
- · Disk brake + V-Brake
- Gear box
- ·Adjustable pedaling position

WITH PEDAL ASSISTANCE




Evolution of the Tribike Hybrid, it maintains the same characteristics of compactness and manageability, adding further eletronic and safety features. Such as reverse gear, the exclusive braking system IBS (Intelligent Braking System) and the ICC system (Intelligent Cruise Control) Freedom and well-being in complete autonomy thanks to pedal assistance.

Notes			



WHEELCHAIR FEATURES			USER'S PERSONAL DATA		
Brand			Name		
Model			Tel.:		
Front tube's diam	neter		E-mail:		
Seat width			Age		Pathology:
Seat depth (E)					
Toe to front tube	distance (A)		DELIVERY ADDRESS		
Ground/knee heig	ght (D)		DLLIVLK	I ADDR	KC33
Length of stretched arms (C)			Delivery to:		
Front tubes angle	(degree)				
Frame length (B)			Postcode	е	
Front tubes width	ı(F)*		Address		
User's weight					
Footrest			Town		
Fixe	d	Internal	County		
Removabl	le	External			

^{*} The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm **All fields are obligatory**



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc.)

Assembly, installation and instruction for use carried out by a Triride specialist technician (on request)



MTW: TREKKING POWER				
Adaptation for MTW Trekking Power Kit purchase (complete and attach dedicated form) no cost				
T0013	Adaptation for MTW Trekking Power(future purchase)	£185.00		
TRIBIKE CONFIGURATION				
Battery				
BS-R-394-S	Battery 36V 9 USB (324Wh)	standard		
BS-R-393-S	Battery 36V 14 USB (504 Wh)	£252.00		
BS-R-395-S	Battery 48V 10.5 USB (504 Wh)	£306.00		
Tyre				
BS-R-054	16" City	standard		
BS-R-051	16" Cross	standard		
BS-R-066	20" City	standard		
BS-R-259	20" Cross	standard		
Cranks				
TR-B-032	120mm	standard		
TR-B-021	140mm	standard		
TR-B-004	155mm	standard		
Handlebar sten	n			
BS-O-221	Fixed 60mm	standard		
TR-B-090	Fixed 80mm 40 degree	standard		
TR-B-091	Fixed 100 mm 40 degree	standard		
BS-C-441	Fixed 120mm 40 degree	standard		
Frame colour				
BS-O-217-N	Black Ral 9005	standard		
BS-O-217-R	Red Ral 3020	standard		
BS-O-217-W	White Ral 9003	standard		
BS-O-217-G	Yellow Ral 1003	£135.00		
BS-O-217-A	Orange Ral 2004	£135.00		
BS-O-217-B	Blue Ral 5017	standard		
BS-O-217	RAL on request	£207.00		



Mudguard		
BS-O-245	Handlescrews - Black	default
BS-O-167	Handlescrews - Blue	£0.00
BS-O-169	Handlescrews - Orange	£0.00
BS-O-166	Handlescrews - Green	£0.00
BS-O-164	Handlescrews - Red	£0.00

OPTIONAL

Attachment K	it for SECOND	WHEELCHAIR - please	e provide d	etails below(prices vary as	per specifics)	
Tilting frame		Connecting arm	s	Clamps		
Brand		Model				
Seat width	cm	Seat depth	cm	Front tubes diameter	mm	

ACCESSORIES				
BS-R-394	Battery 36V 9A USB (324Wh)	£567.00		
BS-R-393	Battery 36V 14 USB (504 Wh)	£702.00		
BS-R-395	Battery 48V 10.5 USB (504 Wh)	£753.00		
BS-O-278	Lateral stand with wheels	£113.00		
BS-O-251	Carry Trailer for fixed frame Rear axle diameter mm Leg	£279.00		
BS-O-214	cover	£19.00		
BS-O-292	Calf support with adaptation to tilting frame	£142.00		
BS-O-114	Arms stopper rings, pair	£16.00		
BS-O-174	Laser customised engraving on frame (on inox and black frame only)	£68.00		
BS-O-156	Digital engraving	£68.00		
BS-O-157	Personalised sticker	£59.00		
BS-O-081	Integrated front light, activated from display	£95.00		
BS-O-232	Light - Front 1200 Lumen	£100.00		
BS-O-198	Light Front USB	£29.00		
BS-O-180	Light - Rear USB	£33.00		
BS-O-015	Rear axle extension for wheelchair (please complete dedicated form)	on request		

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website web www.triride.uk.com

N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product