

Triride®

by Gianni Conte

A NEW FREEDOM



mod. TRIBIKE

Medical Devices Register of the Health Ministry. TUV Certification No. 15-PEL-0010022-TIC

CE



16" / 20"



Manual



9,5 Kg

TECHNICAL FEATURES

EN

Fully adjustable pedal position

Double brake (Disk brake + V-brake)

Aluminium 16" or 20" driving wheel

Tyre: City or Cross

Steel frame

8-speed Shimano Nexus Gear

High resistance painting

Customizable color frame



¹ Weight refers to device excluding arms for wheelchair connection.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.