



TRIBIKE MANUAL

£2,590.00

ORDER ☐

QUOTE ☐

ORDER REFERENCE

DATE

Please complete and send this form by e-mail to: orders@triride.uk

TECHNICAL FEATURES

- 1 - 20" aluminium rim - double chamber
- 2 - Weight: 8 Kg approx.
- 3 - Disc brake + V-BRAKE brake
- 4 - Steel / stainless steel frame
- 5 - Shimano Nexus 8-speed hub
- 6 - Adjustable pedaling position - **N.B.: after each adjustment check the chain length**

USER'S PERSONAL DATA*

Name

Surname

Phone

E-mail

Weight

Age

Pathology

DELIVERY ADDRESS

Delivery to

Street

Town

Postcode

WHEELCHAIR FEATURES

Brand

Model

Front tubes diameter mm

Seat width cm

Seat depth (b) cm

Toe to front tube distance (a)* cm


Ground/knee height (d) cm

Length of stretched arms from the back(c) cm

Footrest ☐ Fixed ☐ Internal ☐

☐ Removable ☐ External ☐

* For wheelchairs with removable footrests, measure as per "a" from toe to front tube.



CONFIGURATION				
Wheel	16" wheel		<input type="checkbox"/>	standard
	20" wheel		<input type="checkbox"/>	standard
Crancks	TR-B-021	140mm (recommended for 16" wheels)	<input type="checkbox"/>	standard
	TR-B-004	155mm (recommended for 20" wheels)	<input type="checkbox"/>	standard
Crown FSA	TR-B-028	42 cogs (recommended Tribike Hybrid)	<input type="checkbox"/>	standard
	TR-B-027	38 cogs (recommended Tribike Manual)	<input type="checkbox"/>	standard
Handlebar Stem	BS-O-221	Fixed 60mm	<input type="checkbox"/>	standard
	TR-B-090	Fixed 80mm 40 degrees	<input type="checkbox"/>	£35.00
	TR-B-091	Fixed 100mm 40 degrees	<input type="checkbox"/>	£35.00
	BS-C-441	Fixed 120mm 40 degrees	<input type="checkbox"/>	£35.00
Tyre	CA-R-120	City	<input type="checkbox"/>	standard
	CA-R-144	Cross	<input type="checkbox"/>	standard
	CA-R-010	Marathon Schwalbe	<input type="checkbox"/>	£14.00
Colour	BS-O-217-N	Black - Ral 9005	<input type="checkbox"/>	standard
	BS-O-217-R	Red - Ral 3020	<input type="checkbox"/>	standard
	BS-O-217-W	White - Ral 9003	<input type="checkbox"/>	£132.00
	BS-O-217-G	Yellow - Ral 1023	<input type="checkbox"/>	£132.00
	BS-O-217-A	Orange - Ral 2004	<input type="checkbox"/>	£132.00
	BS-O-217-B	Blue - Ral 5015	<input type="checkbox"/>	£132.00
	BS-O-217	RAL on request RAL no. _____	<input type="checkbox"/>	£177.00
OPTIONAL				
T-Stand with wheels	BS-O-201		<input type="checkbox"/>	£109.00
Leg cover	BS-O-214		<input type="checkbox"/>	£15.00
Laser customized engraving on the frame	BS-O-174		<input type="checkbox"/>	£68.00
Light - Front 1200 Lumen	BS-O-232		<input type="checkbox"/>	£100.00
Light - Front+USB	BS-O-198		<input type="checkbox"/>	£27.00
Light - Rear+USB	BS-O-180		<input type="checkbox"/>	£32.00
***Prices are intended for the public and before tax				
N. B.:The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product				