

Ventus

The Ventus is a rigid frame wheelchair for active users who demand a lot from their wheelchair. The ventus is a wheelchair with manufactured with sophisticated technology, featuring excellent operating performance and reliability.

Adjusting the centre of gravity, front and rear seat height and back angle is straightforward with the ventus. With quick release rear wheels and a back that folds forward the ventus is easily transferred into a vehicle.

The timeless and efficient design of this agile wheelchair results in an all-round high quality product for everyday use is available in 13 frame colours.

What is special about me?

- Open modern design frame
- Outstanding quality with low rolling resistance enabling wheelchair drivers to use for longer periods before feeling tired or fatigued



Ventus

Standard Model

- Aluminium tube frame with open design
- Folding back rest
- Angle adjustable backrest
- Aluminium side panels
- Wheel camber selectable in 3° increments
- Hollow rims with high pressure tyres
- Knee lever wheel lock
- Rigid axle
- Height and angle adjustable footplate

Options

- Various side panels
- Wheel lock lever extension
- Wheelbase extension
- Selection of push handles
- Straight back tube, fixed
- Various side panels
- Shock absorber
- Frame Pad
- Anti-tipper
- Tip-assist
- Selection of footplates
- Lap Belt
- Various back upholsteries



Height adjustable push handles



Single panel footplate



Swing away anti tipper



Tension adjustable back

Technical Data

Seat width	280-500mm
Seat depth	300-500mm
Front seat height	400-520mm
Rear seat height	330-520mm
Back height	300-500mm
Lower leg length	200-500mm
Width	450-790mm
Length	755-970mm
Height	630-1050mm
Max user weight	100kg (up to SW 440mm)
Transport weights	
Frame	Approx. 8.1kg
Rear Wheel	Approx. 1.9kg
Front frame angle	75°, 85°