



### Tyre Pressures

<b>Model</b>	<b>Front (psi)</b>	<b>Rear (psi)</b>
Buddy/Superlight	36	36
Breeze S4	25	40
Breeze S4 GT	25	30
Breeze S3	40	40
Mystere	50	50
Breeze Midi	25	25
Sonet	40	40
Supersport	35	40
Vita 4/Sport	32	32
Vita Lite	30	30
Vita Midi	30	30
Vita X	25	25
Minimo Plus	20	20
Maximo	30	30