



WOLTURNUS

# W5 Junior

Tailor-made rigid wheelchair  
for children



# Dynamic, stable and light

## Ideal for the active child

W5 Junior is a customized active rigid wheelchair, designed for optimum synergy between strength and weight.

The W5 Junior's backrest is fully welded and it has an adjustable rear axle, which allows the continuous adaption of the balance point. However, the wheelchair is also available with a welded rear axle. Then, the wheelchair has no moving parts and is therefore extremely rigid and robust, which means that all of the child's energy is harnessed for propulsion.

The ergonomic seat has an adjustable strap seat as well as a strap back.

This ensures an ergonomic, upright posture, optimum pressure relief and freedom of movement. The seat depth can be adjusted continuously. With the help of an extension-kit, it is possible to extend the seat width by up to 3 cm, as your child grows.

The castors are positioned out in front of the foot plate, which ensures high stability and low rolling resistance.

There are many accessories, such as push handles or spokeguards, available to choose from.

The W5 Junior can be ordered in different colors.

### FEATURES

- › Customized rigid wheelchair
- › Ideal for the active child
- › Extremely sturdy
- › Fully welded or adjustable versions available
- › Very easy to manoeuvre
- › Multi-adjustable footrest
- › High stability and low weight
- › Crash-tested and approved for fastening in a car
- › Propulsion weight from 5 kg
- › User weight max. 70 kg



#### › REAR AXLE

You can choose between a fixed rear axle and an adjustable rear axle.



#### › ERGONOMIC SEAT

Ergonomic seat with adjustable strap back and strap seat.



#### › CASTORS

The position of the castors provides extra stability and a low rolling resistance.