

<p>Base Price - £2846.00</p> <p>Valid from 01.01.24 to 31.12.24 Please ensure you are using the latest order form</p>	<p>Order no.:</p> <p>-----</p>
--	--------------------------------

<input type="checkbox"/> Provide estimate only	<input type="checkbox"/> Order
--	--------------------------------

Billing address/customer number

Delivery address

Sales person	<p>Name</p> <p>.....</p> <p>Date Signature</p>
--------------	---

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Wingus

Power wheelchair

Max. load capacity: 120 kg
Article number: 490E163=00011



Standard equipment

- nVR2 control unit (40 A)
- Max. speed: [6 km/h or 7,2 km/h]
- Standard control holder, mounting side right
- Batteries: [AGM or Gel]
- Max. distance range 25 km
- Battery charger: [5 A or 10 A]
- Caster wheels 8" (grey)
- Drive wheels 12,5" (grey)
- Puncture-proof wheels
- Seat incl. seat cushion (W x D x H: 45 x 51 x 5 cm)
- Back angle adjustment mechanical (0°/10°/20°)
- Flip-up arm support
- Height-adjustable leg support with swing-up foot plate, angle adjustable (-5°/0°/+5°)
- LLL adjustable: 38 - 48 cm
- Lap belt with clip closure
- Frame colour: Black (CDP)
- Empty weight: 57 kg

Article number		Article	Quantity	Price
490E163=00011	<input type="checkbox"/>	Wingus		£2846.00