



WOLTURNUS

Dalton

Rigid wheelchair with
multifunctional footrest



Functional and versatile

The Dalton is a rigid wheelchair for active users, who want a multi-adjustable and detachable footrest.

The Dalton is a rigid wheelchair with separate flip-up foot plates, which are adjustable in angle and height. It has swing-away, detachable leg supports which can be mounted or removed easily. This allows an easy and save transfer into and out of the wheelchair.

If it should be necessary to elevate the legs, it is possible to equip the Dalton with an elevating footrest on just one or both sides.

The Dalton is customized to suit the individual user's level of activity, size and strength.

The traditional model can be used up to a max. user weight of 160 kg. For strongly-built users, the Dalton can be designed as an extra strong version, which allows a max. user weight of 250 kg. For users with max. 80 kg of body weight, Wolturnus offers a lighter alternative to the traditional Dalton version. The Dalton Light has slimmer frame tubes, leading to a more elegant look.

By choosing a steplessly adjustable balance point and an angle and height-adjustable folding back, you can continuously adapt the seating posture. To reduce the weight, the Dalton can also be ordered fully welded.

Finally, the Dalton is also available with a particularly low seat height. The Dalton L is ideal for users who want to 'walk the chair forward'.


FEATURES

- › Separate, flip-up foot plates that are angle and height-adjustable
- › Swing-away, detachable leg supports
- › Elevating footrest
- › Detachable or fixed side panels
- › Angle and height-adjustable folding back (optional)
- › Available with Wing Back ILSA
- › Available with low seat height (Dalton L)
- › Frame made of aluminium 7020
- › Weight from 10 kg
- › User weight max. 160 kg (250 kg)



› FOOTREST

Separate, flip-up foot plates, which are angle and height-adjustable. The swing-away function allows an easy transfer.



› ELEVATING FOOTREST

You can choose an elevating footrest on one or both sides.



› DALTON LIGHT

The Dalton Light comes with slimmer frame tubes, resulting a more elegant look. Max. user weight 80 kg.