



WOLTURNUS

# HAWK

Live Your Life Actively







# HAWK

Your Life. Your Choices.

## THE HAWK

The hawk is a fast middle-sized bird of prey that rarely stays idle for longer periods of time. It is a big powerful bird with short, wide wings and a long tail. The hawk is often seen flying low above the ground and gliding high up in the air.







# HAWK

HAWK is created for the ultimate unfoldment of your life's aspects. A compact, sturdy and unique wheelchair for active use with a solid frame, an impressive drivability and a vast number of customization options. It is adjustable to fit any wheelchair user's dreams and challenges.

# HAWK IS CREATED FOR YOU

## Your life. Your choices. Your needs

The HAWK is built to fit you, for the life you live, for your choice of accessories, your unique seating posture and your preferences. The HAWK is elegant and handy, which is capable of dealing with the challenges of every-day-life's challenges, indoors and outdoors. The HAWK adjusts intuitively - it is modifiable to fit your lifestyle.

Your HAWK is unique since you can customize it with the colors and components reflecting your choices. This wheelchair is easy to order and it has short delivery times - however the HAWK is anything but a standard item from the warehouse shelf helping you to overcome your mobility handicap. On contrary, it is your faithful companion throughout the unfoldment of your life.



### FACTS AND FIGURES:

- › Active rigid wheelchair
- › 7020 aluminium box-frame
- › Power-coated frame
- › Ergonomic seat
- › Colorful components
- › Adjustable strap seat and strap back
- › Adjustable seat height and seat depth
- › Folding back with adjustable back angle
- › Weight from 7kg
- › Maximum slope, using parking brakes: 7°
- › Maximum weight of the user: 120 kg
- › Easy to fit into a car
- › CE-approved & TÜV-certified



ALL INFORMATION ABOUT THE HAWK ARE ALSO AVAILABLE ON  
**WOLTURNUS.COM**



The best thing about my  
Wolturnus wheelchair is that  
I sometimes forget it exists.





# A WHEELCHAIR THAT ALWAYS FITS

The HAWK is made for all active wheelchair users who need an adjustable wheelchair for daily use. The wheelchair is very precisely customized to fit your measures - actually, down to the last millimeter.

The HAWK has features such as: an adjustable balance point, modifiable seat height and depth, adjustable strap seat and adjustable strap back. Furthermore, it has an adjustable back angle and adaptable back height.



# UNIQUE DESIGN



## » BOX FRAME

A closed box-frame made of 7020 aluminium unites optimal strength and stability with a modern and elegant design.



## » DIBOND SIDEGUARDS

Detachable Dibond sideguards, which can be removed.



## » FOLDING BACK

A specially designed folding back, which locks when folded. This makes it easy to get the chair in and out of a vehicle. The folding back is adjustable by angle.



## » ERGONOMIC SEAT

The rear part of the seat is flat and its front inclined to an upright position and minimizes the risk of pressure.



The HAWK has unique features, united in an elegant, practical and sturdy design that suit your specific needs.



oved within seconds



► **ADJUSTABLE FOOTPLATE**

An ultralight and elegant footplate, with adjustable height and adjustable angle, as well as anti-slip coating.



es. The seat stabilizes the pelvis towards an  
e ulcers.



► **PUSH WHEEL LOCKS**

As a standard, the HAWK is equipped with reliable push-wheel locks.

# YOUR HAWK 100% YOUR CHOICES

## Create your unique HAWK with components of your choice

These are only a few out of a huge amount of options you have, creating your individual wheelchair

Go to the last pages and read all about the standard configurations as well as additional accessories



### » COLORED COMPONENTS

Choose between 7 different component colors



### » FIBER PLASTIC WHEEL LOCKS


Ultralight scissor wheel locks, for avoiding additional weight



### » WING BACK ILSA

A specially designed, patent pending, multi-adjustable back rest system for individual modification of the sitting position





to "spice up" your front forks, sideguard poles, back brackets and footrest tubes



› **SPINERGY REAR WHEELS**

Spinergy wheels with spokes in your favorite color



› **CURVE PUSH RIM**

Special pushrims, with ergonomically optimized profile



› **ADJUSTABLE PUSH HANDLES**

Heightadjustable push handles with rubber coating



› **CARBON SIDEGUARDS WITH MUDGUARDS**

Detachable Carbon sideguards with mudguards.



# JUST IN TIME

The HAWK is easy to use - and easy to order. Although you create your HAWK based on your personal taste and needs, it can be delivered very fast.

## 5 easy steps to order your HAWK

- 1** Contact an approved Wolturnus partner
- 2** Get help in choosing your seat width and other important measures
- 3** Choose your favorite accessories
- 4** Fill in a order form together with your Wolturnus partner
- 5** Your specific HAWK is delivered in a matter of days through your Wolturnus partner





## FOR A LIFE ON THE MOVE

The HAWK is easy to disassemble and transport in a car. The folding back locks when it is folded together. This makes it easy to get the wheelchair into a car while you grab the handle on the back.



Technical Data

Measurements	SB 350	SB 375	SB 400	SB 425	SB 450	SB 475	SB 500
Seat width (SW) (mm)	350	375	400	425	450	475	500
Total width (TW) (mm)	515	540	565	590	615	640	665
Seat height front (mm)	460-520	460-520	460-520	460-520	460-520	460-520	460-520
Seat height rear (mm)	390-480	390-480	390-480	390-480	390-480	390-480	390-480
Back height (BH) (mm)	300-500	300-500	300-500	300-500	300-500	300-500	300-500



ALL INFORMATION ABOUT THE HAWK ARE ALSO AVAILABLE ON  
**WOLTURNUS.COM**



## Technical Data

Frame	Standard Configuration	Alternatives/ Additional Accessories
<b>Total width</b>	Seat width +19 cm , at 1° camber angle	
<b>Seat width</b>	35, 37.5, 40, 42.5, 45, 47.5 or 50 cm	
<b>Seat depth</b>	35-41 / 39-45 / 43-50	
<b>Balance point</b>	Steplessly adjustable, between 2-9 cm	
<b>Seat height front/back</b>	48/42 cm	adjustable from 46-52 cm / 39-48 cm
<b>Front frame</b>	UV-shaped	
<b>Lower leg length</b>	adjustable 39 - 45 cm (std.)	adjustable 36 - 50 cm
<b>Footplate</b>	adjustable by angle, 15 cm deep with anti-slip coating	
<b>Footrest</b>	steplessly height adjustable	
<b>Push handles</b>	None	Height adjustable or foldable push handles
<b>Wheel locks</b>	Push wheel locks	Pull-activated wheel locks, horizontal push wheel locks, scissor wheel locks, scissor wheel locks (fiber), handle extension
<b>Anti-tipper</b>	None	Wing anti-tipper = swing-away option (left or right side)
<b>Frame colors</b>	Powder coating: white, black or grey	
<b>Component colors</b>	Purple, green, red, blue, gold, silver, black	
<b>Colored components</b>	Front forks, sideguard pole, back brackets, footrest tubes	

Wheels	Standard Configuration	Alternatives/ Additional Accessories
<b>Rear Wheels</b>	24" standard cross-spoke wheels with quick-release axles	22" / 25" : standrad cross spoke wheels 22" / 25" Spinerger (12/18 spokes - black/white/blue/yellow/red/green/orange/pink)
<b>Pushrims</b>	24" aluminium silver-anodized	22" / 25" aluminium silver-anodized 22" / 24" / 25" : Titanium, Friction, Surge LT, plastic, steel, Maxgrepp Ergo Para, Curve L TG, B-Braver
<b>Tyres</b>	24" Standard tyres (grey)	22" / 25" Standard tyres (grey), 24" Puncture resistant (massive), 22" / 24" / 25" Schwalbe Marathon Plus (black) 24" / 25" Schwalbe Marathon Plus (grey)
<b>Camber Angle</b>	3°	1°
<b>Front Wheels</b>	4"/34 mm fiberplast castors (black, blue, red, white, silver)	5" / 6" fiberplast castors (black) 4" / 5" : Frog Legs (Aluminium) in various colors, 6" black

Positioning	Standard Configuration	Alternatives/ Additional Accessories
<b>Back Frame</b>	Foldable with back lock, height adjustable 35-40 cm (std.), adjustable by angle ±4°	height adjustable between 30-50 cm Wing Back ILSA (3,5, 7 cm depth)
<b>Seating Frame</b>	Ergonomic seat	
<b>Adjustable strap seat</b>	Black nylon upholstery with adjustable seating depth	
<b>Adjustable strap back</b>	Black, double-layered upholstery with foam filling, breathable and ventilating	
<b>Sideguards and armrest</b>	Foldable or detachable Dibond sideguards	Carbon without mudguards, (detachable/ foldable), Carbon with mudguards (detachable), ABS sideguards with armrest (detachable), sideguards with thermal protection



**WOLTURNUS**

Wolturnus A/S | Skalhuse 31 | DK-9240 Nibe | T. +45 96 71 71 70 | E. [info@wolturnus.dk](mailto:info@wolturnus.dk) | [www.wolturnus.com](http://www.wolturnus.com)

PARTNER:

**ottobock.**



Recare  
13 Lodge Road  
Hanborough Business Park  
Long Hanborough  
Oxfordshire  
OX29 8LJ

[www.recare.co.uk](http://www.recare.co.uk)  
Tel. 01993 880200